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SUMMARY AND EXPLANATION

PYR reagent is reagent used for generation of colour in test PYR in diagnostic kits. Detailed description of work with reagent PYR and interpretation of result of test is part of particular kit manual.

PRINCIPLE

Substrate L – pyrrolidonyl – b – naftylamid is transformed by acting of L – pyrrolidonylmarylamidase to colourless L – pyrrolidon and b - naftylamin. Free b – naftylamin is visualized by PyR reagent by creation of red coloured Shiffs base.

DIAGNOSTIC KIT CONTAINS

- PYR A (Ref. 3003)
- PYR reagent (Ref. 3003)
- dropper
- information leaflet

PREPARATION OF REAGENT

Pour content of glass bottle to plastic bottle, and then add dropper and lid to bottle. PYR reagent is ready for use now.

WARNINGS AND SAFETY PRECAUTIONS

- **For *in vitro* diagnostics use.**
- **For professional use only**
- Follow the instructions exactly!
- All samples and inoculated products have to be considered as potentially infectious and obvious safety regulations of your country must be respected while handling.
- Do not use overexpired product.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of overexpired or devaluated reagent liquidate according to regulations of your country. Packaging materials are recyclable.

R 36/R38 Irritating to eyes and skin.

S 37/39 Wear suitable gloves and eye/face protection.

S 26 In case of contact with eyes rinse immediately with plenty of water and seek medical advice.

S 45 In case of accident or if you feel unwell seek medical advice immediately (show the label when possible).

STORAGE

PYR reagent can be stored in dark in temperatures between +2 and +8°C. Expiration date is showed on label of each package.

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