A Commitment to Diagnostics

# CAST<sup>®</sup> Allergens

for

CAST<sup>®</sup> ELISA  
Flow CAST<sup>®</sup>  
Flow CAST<sup>®</sup> highsens

Drug Allergens

<b>Introduction</b>	<b>4</b>	<b>Antiseptics</b>	<b>27</b>
<b>Literature</b>	<b>14</b>	BAG2-CCHX Chlorhexidine	27
<b>Adjuvants</b>	<b>15</b>	<b>Contrast Media</b>	<b>28</b>
BAG2-CE466 Carboxymethylcellulose	15	BAG2-CIOBI Iobitridol	28
BAG2-CPEG PEG-4000	15	BAG2-CIODI Iodixanol	28
<b>Analgesics</b>	<b>16</b>	BAG2-CIOHE Iohexol	28
BAG2-C51 Lys-Aspirin	16	BAG2-CIOME Iomeprol	29
BAG2-C52 Diclofenac	16	BAG2-CIOPA Iopamidol	29
BAG2-C53 Ibuprofen	16	BAG2-CIOPR Iopromide	29
BAG2-C54 Indomethacin	17	BAG2-CIOXA Ioxaglate	30
BAG2-C54F Indomethacin	17	BAG2-CIOXI Ioxitalamate	30
BAG2-C55 Acetaminophen	17	<b>Local Anesthetics</b>	<b>31</b>
BAG2-C56 Mefenamic Acid	18	BAG2-CART Articaine	31
BAG2-C57 Phenylbutazone	18	BAG2-CBUP Bupivacaine	31
BAG2-C58 Propyphenazone	18	BAG2-CMEP Mepivacaine	31
BAG2-C59 Dipyron / Metamizole	19	BAG2-CLID Lidocaine	32
BAG2-CNAP Naproxen	19	<b>Myorelaxants</b>	<b>33</b>
BAG2-C114 Salicylic acid	19	BAG2-CATR Atracurium	33
<b>Antibiotics</b>	<b>20</b>	BAG2-CCAT Cis-Atracurium	33
BAG2-C1 Penicillin G	20	BAG2-CMIV Mivacurium	33
BAG2-C11 Benzylpenicilloyl-Polylysine (PPL)	20	BAG2-CPAN Pancuronium	34
BAG2-C12 Minor Determinant Mix (MDM)	20	BAG2-CPRO Propofol	34
BAG2-C2 Penicillin V	21	BAG2-CROC Rocuronium	34
BAG2-C3 Cephalosporin C	21	BAG2-CSUX Suxamethonium	35
BAG2-C31 Cefamandole	21	BAG2-CVEC Vecuronium	35
BAG2-C32 Cefazoline	22	BAG2-CNBX Neuromuscular Blockers Mix	36
BAG2-C33 Cefuroxime	22	<b>Proton Pump Inhibitors</b>	<b>37</b>
BAG2-C34 Cefaclor	22	BAG2-COME Omeprazole	37
BAG2-C35 Ceftriaxone	23	BAG2-CESO Esomeprazole	37
BAG2-C36 Levofloxacin	23	BAG2-CPANT Pantoprazole	37
BAG2-C61 Sulfamethoxazole	23	BAG2-CLAN Lansoprazole	38
BAG2-C62 Trimethoprim	24	BAG2-CRAB Rabeprazole	38
BAG2-C75 Tetracycline	24		
BAG2-C81 Ciprofloxacin	24		
BAG2-C203 Ampicillin	25		
BAG2-C204 Amoxicillin	25		
BAG2-CCLA Clarithromycin	25		
BAG2-CCLAV Clavulanic acid	26		
BAG2-CRIF Rifampicin	26		

## Intended Use

BÜHLMANN CAST® Allergens are optimized for *in vitro* Basophil Activation using BÜHLMANN CAST® assays: CAST® ELISA (EK-CAST), Flow CAST® (FK-CCR) and Flow CAST® highsens (FK-HSBAT).

CAST® Allergens are not included in BÜHLMANN CAST® assays. They are provided upon request. Refer to the allergen list.

## Instruction for Use

The following instructions are intended for the following groups of drug allergens and mixes of them (refer to TOC page 3).

- Adjuvants
- Analgesics
- Antibiotics
- Antiseptics
- Contrast Media
- Local Anesthetics
- Myorelaxants
- Proton Pump Inhibitors (PPI)

CAST® drug allergens are delivered lyophilized. Before usage in CAST® assays, they are to be reconstituted with test specific Stimulation Buffer: CAST® ELISA: B-CAST-STB, Flow CAST®: B-CCR-STB and Flow CAST® highsens: B-CCR-CSB.

- Add 250 µl Stimulation Buffer into the vial and vortex until the lyophilized material has completely dissolved.
- For some allergens we recommend an additional dilution with Stimulation Buffer (refer to the paragraph „Additional dilutions“ in the allergen specifications).
- E.g. an „additional dilution“ of 1:5 is carried out by adding 160 µl Stimulation Buffer to 40 µl allergen solution.
- For stimulation of cells, use these dilutions according to the respective CAST® IFU.

Remark: Some individuals react positive either at high or at low allergen concentrations, whereas others react positive over a broad range of concentration.

## Storage and Stability

Unopened CAST® Allergens are to be stored at the temperature specified on the label until to the expiration date.

**IMPORTANT:** Reconstituted or diluted allergens must not be stored and used again!

## Optimal time period for stimulation tests

For optimal results, CAST® assays should be run between 3 and 12 weeks after an allergic reaction to the presumed allergen [Lit.8], although specific basophil reaction is often being persistent over a longer time period.

Blood samples must be drawn for CAST® assays before carrying out a skin test or *in vivo* provocation with the presumed allergen as *in vivo* exposure to an allergen can cause an activation of the leucocytes.

## Limitations

A negative CAST® assay result for a specific allergen does not exclude potential (even serious) clinical symptoms. Individuals who developed „allergic“ reactions towards protein or drug allergens in the past and whose CAST® assay was negative, should thus be verified with additional methods like the *in vivo* provocation or skin test before a drug or allergen is administered.

## Concentrations and cut-off

Concentrations and cut-off of these allergen groups were optimized regarding maximum specificity und sensitivity. Concentrations and cut-offs of individual allergens are reported in the specifications of the allergens (refer to TOC page 3).

For Flow CAST®, BÜHLMANN recommends the application of the stimulation index (SI) as an additional diagnostic criterion. The SI is defined as the ratio of allergen specific basophil activation and background activation.

## Remarks:

- In order to apply allergen specific cut-offs, QC criteria published in the CAST® instructions must be fulfilled.
- Cut-offs should serve only as recommendations. Clinically defined cut-offs should be established by each laboratory or through clinical studies.

## Recommendations how to use CAST® assays in the diagnosis of drug allergies

The diagnosis of drug allergies is a complex matter including heterogeneous symptoms and a diverse pathophysiology. Thus the following recommendations should be considered when using CAST® assays in order to diagnose drug allergies. These recommendations should be considered only as general guidelines.

---

### Clinical pictures, justifying the application of CAST® assays

Direct "allergic" reactions *in vivo* after exposure to drugs or chemicals can either be IgE or non-IgE mediated (pseudo-allergic reactions) stimulation of effector cells (e.g. basophils, mast cells, and eosinophils). CAST® assays may serve as a model for the above mentioned reactions by which the individuals' leucocytes are stimulated *in vitro*. As for CAST® assays pure allergens and a defined but artificial buffer system are applied, the method can only serve as a simplified model of the real *in vivo* conditions.

CAST® assays were optimized for the detection of immediate type allergic reactions (Type I) and pseudo allergic reactions. CAST® assays are not applicable for the diagnosis of delayed type allergies (e.g. Type IV reactions). Thus, these assays should be applied only for immediate type allergic reactions, especially if one of the following symptoms is detected:

- Anaphylactic or anaphylactoide reactions
- Rhino conjunctivitis
- Asthma bronchiale
- Urticaria / Angioedema
- Gastrointestinal reactions

### Clinical pictures that do not justify the application of CAST® assays

CAST® assays should not be used for the diagnosis of T-cell mediated reactions. Based on the pathogenetical mechanism it is not recommended using CAST® assays routinely for the following patient groups:

- Maculopapular and pustular exanthema
- Vasculitis
- Allergic contact dermatitis
- Bullous exanthemas such as Stevens-Johnsons Syndrome and Lyell's disease

### Use of allergens not available at BÜHLMANN

Allergens which are not available at BÜHLMANN can be used in CAST® assays. The following restrictions should be observed:

- No allergens bound to matrices (solid or liquid)
- No allergen preparations containing sulfido-leukotrienes
- No allergen preparations containing cytotoxic ingredients (stabilizers, preservatives), like e.g. phenol, glycerol, sodium azide or merthiolate (thimerosale)

An instruction how to make allergen extracts which are not available at BÜHLMANN, is available upon request. Ask your local distributor or contact BÜHLMANN Laboratories AG directly.

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1. Eberlein B, León Suárez I, Darsow U, Ruëff F, Behrendt H, Ring J: A new basophil activation test using CD63 and CCR3 in allergy to antibiotics. *Clin & Exp Allergy* 2010, 40: 411-8.
2. Garcia-Ortega P, Marin A, Quilez E, Rodriguez-Leor O: New tools for allergy diagnosis and clinical advice after immediate reactions to ICM. *Br J Radiol* 2010, 83: 632-33.
3. De Weck AL, Sanz ML, Gamboa PM, Jermann JM, Kowalski M, Medrala W, Sainte-Laudy J, Schneider MS, Weber JM, Wolanczyk-Medrala A: Nonsteroidal anti-inflammatory drug hypersensitivity syndrome. A multicenter study II. Basophil activation by nonsteroidal anti-inflammatory drugs and its impact on pathogenesis. *J Investig Allergol Clin Immunol* 2010, 20: 39-57.
4. Hausmann OV, Gentinetta T, Bridts CH, Ebo DG: The basophil activation test in immediate-type drug allergy. *Immunol Allergy Clin N Am* 2009, 29: 555-66.
5. Dewachter P, Nicaise-Roland P, Kalaboka S, Lefèvre J, Chollet-Martin S: Anaphylaxis to amidotrizoate proved by skin testing and flow cytometry-based basophil activation test. *Allergy* 2009, 64: 501-2.
6. De Weck AL, Sanz ML, Gamboa PM, Aberer W, Sturm G, Bilo MB, Montroni M, Blanca M, Torres MJ, Mayorga L, Campi P, Manfredi M, Drouet M, Staint-Laudy J, Romano A, Merk H, Weber JM, Jermann TM, members of ENDA (European Network for Drug Allergy): Diagnosis of immediate-type  $\beta$ -Lactam allergy in vitro by flow-cytometric basophil activation test and sulfidoleukotriene production: a multicenter study. *J Investig Allergol Clin Immunol* 2009, 19: 91-109.
7. Steiner U, Gentinetta T, Hausmann O: IgE-mediated anaphylaxis to intraarticularn Glucocorticoid preparations. *AJR Am J Roetgenol* 2009, 193:156-7.
8. Gómez E, Blanca-Lopez N, Torresz MJ, Requenaz G, Rondonz C, Canto G, Blanca M, Mayorga C: Immunoglobulin E-mediated immediate allergic reactions to dipyrone: value of basophil activation test in the identification of patients. *Clinical & Experimental Allergy* 2009, 39: 1217-1224.
9. Kvedariene V, Kamey S, Ryckwaert Y, Rongier M, Bousquet J, Demoly P, Arnoux B: Diagnosis of neuromuscular blocking agent hypersensitivity reactions using cytofluorimetric analysis of basophils. *Allergy* 2006; 61: 311-5.
10. Scherer K, Studer W, Figueiredo V, Bircher A: Anaphylaxis to isosulfan blue and cross-reactivity to patent blue V: Case report and review of the nomenclature of vital blue dyes. *Ann Allergy Asthma Immunol* 2006, 96: 497-0.
11. Erdmann SM, Ventocilla S, Moll-Slodowy S, Sauer I, Merk HF: Basophilenaktivierungstests in der Diagnostik von Arzneimittelreaktionen. *Hautarzt* 2005, 56: 38-43.
12. Sanz ML, Gamboa P, de Weck AL: A new combined test with flowcytometric basophil activation and determination of sulfidoleukotrienes is useful for in vitro diagnosis of Hypersensitivity to Aspirin and other nonsteroidal anti-inflammatory drugs. *Int Arch Allergy Immunol* 2005, 136: 58-72.
13. Sudherr PS, Haal JE, Read GF, Rowbottom AW, Williams PE: Flow cytometric investigation of peri-anaesthetic anaphylaxis using CD63 and CD203c. *Anaesthesia* 2005, 60: 251-56.
14. Gamboa P, Sanz ML, Caballero MR Urrutia I, Antepara I, Esparza R, de Weck AL: The flowcytometric determination of basophil activation induced aspirin and other non-steroidal anti-inflammatory drugs (NSAIDs) is useful for in vitro diagnosis of the NSAID hypersensitivity syndrome. *Clin Exp Allergy* 2004, 34: 1448-57.
15. Torres MJ, Padial A, Mayorga C, Fernandez T, Sanchez-Sabate E, Cornejo-Garcia JA, Antunez C, Blanca M: The diagnostic interpretation of basophil activation test in immediate allergic reactions to betalactams. *Clin Exp Allergy* 2004, 34: 1768-75.
16. Gamboa PM, Sanz ML, Caballero MR, Antépara I, Urrutia I, Jáuregui I, González G, Diéguez I, De Weck AL: Use of CD63 expression as a marker of in vitro basophil activation and leukotrienes determination in metamizole allergic patients. *Allergy* 2003, 58: 312-7.
17. Sanz ML, Gamboa PM, De Weck AL: Clinical evaluation of in vitro tests in Diagnosis of immediate allergic reactions to  $\beta$ -lactam antibiotics. *ACI International* 2002, 14: 185-93.
18. Monneret G, Benoit Y, Debard AL, Gutowski MC, Topenot I, Bienvenu J: Monitoring of basophil activation using CD63 and CCR3 in allergy to muscle relaxant drugs. *Clin Immunol* 2002, 102: 192-9.
19. Sanz ML, Gamboa PM, Antepara I, Uasuf C, Vila L, Garcia-Aviles C, Chazot M, De Weck AL: Flow cytometric basophil activation test by detection of CD63 expression in patients with immediate-type reactions to betalactam antibiotics. *Clin Exp Allergy* 2002, 32: 277-86.
20. Monneret G, Benoit Y, Gutowski MC, Bienvenu J: Detection of basophil activation by flow cytometry in patients with allergy to muscle-relaxant drugs. *Anesthesiology* 2000, 92: 275-7.
21. Abuaf N, Rajoely B, Ghazouani E, Levy DA, Pecquet C, Chabane H, Leynadier F: Validation of a flow cytometric assay detecting in vitro basophil activation for the diagnosis of muscle relaxant allergy. *J Allergy Clin Immunol* 1999, 104: 411-8.

Order code	Name	
<b>BAG2-CE466</b>	<b>Carboxymethylcellulose</b>	sodium salt

Trade Names Cellulose glycolate; Cethylose; CMC; Carmethose; Cel-0-Brandt; Glykocellon; Carbose D; Xylo-Mucine; Tylose MGA; Cellolax; Polycell

Content per vial	5 µg	Conc after reconst.	20 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	4 µg/ml	4.5 µg/ml	4.5 µg/ml
Cut-off	≥150 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

Order code	Name	
<b>BAG2-CPEG</b>	<b>PEG-4000</b>	H(OCH <sub>2</sub> CH <sub>2</sub> ) <sub>n</sub> OH MW 3500 - 4500

Trade Names Macrogol 4000; Polyethylenglycol MW 4000

Content per vial	1 mg	Conc after reconst.	4 mg/ml
Additional dilutions	1:10		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	800 µg/ml	910 µg/ml	910 µg/ml
Cut-off	≥ 50 pg/ml	≥15% CD63; SI*≥2	≥15% CD63+CD203c

<b>Order code</b>	<b>Name</b>	$C_{15}H_{22}N_2O_6$
<b>BAG2-C51</b>	<b>Lys-Aspirin</b>	MW 326.4

Trade Names	Acetophen; Acetosal; Acylpyrin; Alcacyl; Asatard; Aspidol; Aspro; Caprin; Colfarit; Contrheuma; Delgesic; Duramax; Ecotrin; Endydol; Entrophen; Rhodine; Salacetin; Salcetogen; Saletin; Solpyron; Venopirin		
Content per vial	625 µg	Conc after reconst.	2.50 mg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	500 µg/ml	570 µg/ml	570 µg/ml
Cut-off	≥90 pg/ml	≥5% CD63; SI*≥2	≥10% CD63+CD203c

<b>Order code</b>	<b>Name</b>	$C_{14}H_{10}Cl_2NNaO_2$
<b>BAG2-C52</b>	<b>Diclofenac</b>	MW 318.1

Trade Names	Allvoran; Benfofen; Cataflam; Diclo-Phlogont; Diclo-Puren; Diclorem; Dolobasan; Duravolten; Ecofenac; Effekton; Novapixina; Primofenac; Prophenatin; Rhumalgen; Valetan; Voldal; Voltaren; Voltarol; Xenid		
Content per vial	6.25 µg	Conc after reconst.	25 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	5 µg/ml	5.7 µg/ml	5.7 µg/ml
Cut-off	≥40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{13}H_{18}O_2$
<b>BAG2-C53</b>	<b>Ibuprofen</b>	MW 206.3

Trade Names	Alges; Algifo; Brufen; Dismenol; Dentigoa; Dolgirid; Dolgit; Dolocyl; Epobron; Femadon; Grefen; Ifufen; Ibumetin; Ibutad; Ibutop; Nurofen; Optifen; Saridon; Sonotryl; Spedifen; Tabalon; Urem		
Content per vial	25 µg	Conc after reconst.	100 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	20 µg/ml	22.7 µg/ml	22.7 µg/ml
Cut-off	≥50 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>		$C_{19}H_{16}ClNO_4$
<b>BAG2-C54</b>	<b>Indomethacin</b>		MW 357.8

Trade Names	Elmetacin; Indocin; Indometacin; Indophtal; Indo-TOP; Indocolir; Indomet; Inflam; Protaxon; Reusin; Rheublamin; Sigadoc; Sjoebo		
Content per vial	2.5 µg	Conc after reconst.	10 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	2.0 µg/ml	2.27 µg/ml	n.a.
Cut-off	≥40 pg/ml	≥5% CD63; SI*≥2	n.a.

<b>Order code</b>	<b>Name</b>		$C_{19}H_{16}ClNO_4$
<b>BAG2-C54F</b>	<b>Indomethacin</b>		MW 357.8

Trade Names	Elmetacin; Indocin; Indometacin; Indophtal; Indo-TOP; Indocolir; Indomet; Inflam; Protaxon; Reusin; Rheublamin; Sigadoc; Sjoebo		
Content per vial	0.5 µg	Conc after reconst.	2 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	n.a.	n.a.	0.45 µg/ml
Cut-off	n.a.	n.a.	≥10% CD63+CD203c

<b>Order code</b>	<b>Name</b>		$C_8H_9NO_2$
<b>BAG2-C55</b>	<b>Acetaminophen</b>		MW 151.2

Trade Names	Acetalgin; Azur; Ben-u-ron; Boxacin; Boxonal; Buscopan; Captin; Cervopyrin; Coffalon; Dafalgan; Dolereduct; Dolevar; Enelfa; Epsimol; Exedrin; Febrifant; Fibrex; Gelonida; Panadol; Perfalgan; Treupel; Zolben		
Content per vial	2.5 µg	Conc after reconst.	10 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	2.0 µg/ml	2.27 µg/ml	2.27 µg/ml
Cut-off	≥60 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{15}H_{15}NO_2$
<b>BAG2-C56</b>	<b>Mefenamic Acid</b>	MW 241.3

Trade Names Mefenacid; Mefenamin; Mephadolol; Ponstan; Spiralgin

Content per vial	6.25 µg	Conc after reconst.	25 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	5 µg/ml	5.7 µg/ml	5.7 µg/ml
Cut-off	≥60 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{19}H_{20}N_2O_2$
<b>BAG2-C57</b>	<b>Phenylbutazone</b>	MW 308.4 sodium salt

Trade Names Ambene; Artrizin; Azolid; Bizolin; Butacote; Butadion; Butapirazol; Butadiona; Butatron; Butoz; Butazolidin; Buzon; Ecobutazone; Equipalazone; Exrheudon N; Fenibutol; Intrabutazone; Intrazone; Mepha-Butazon; Phenyzene

Content per vial	25 µg	Conc after reconst.	100 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	20 µg/ml	22.7 µg/ml	22.7 µg/ml
Cut-off	≥80 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{14}H_{18}N_2O$
<b>BAG2-C58</b>	<b>Propyphenazone</b>	MW 230.3

Trade Names Amtalidin; Budiro; Causyth; Cibalgina; Cistalgin; Eufibron; Isoprochin; Optalidon Saridon Spasmo

Content per vial	12.5 µg	Conc after reconst.	50 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	10 µg/ml	11.5 µg/ml	11.5 µg/ml
Cut-off	≥40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

Order code	Name	C <sub>13</sub> H <sub>16</sub> N <sub>3</sub> NaO <sub>4</sub> S · H <sub>2</sub> O	
<b>BAG2-C59</b>	<b>Dipyrone / Metamizole</b>	MW 351.4	
		sodium salt	
Trade Names	Analgina; Analgin; Berlosin; Inalgon; Metalgin; Minalgin; Nopain; Novalgin; Novalgina; Novaminsulfon; Nolotil; Sulfonovin; Vetalgin		
Content per vial	12.5 µg	Conc after reconst.	50 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	10 µg/ml	11.5 µg/ml	11.5 µg/ml
Cut-off	≥50 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

Order code	Name	C <sub>14</sub> H <sub>13</sub> NaO <sub>3</sub>	
<b>BAG2-CNAP</b>	<b>Naproxen</b>	MW 252.2	
		sodium salt	
Trade Names	Alacetan; Aleve; Anaprox; Apranax; Dolormin; Dysmenalgit; Mobilat; Naprobene; Proxen		
Content per vial	100 µg	Conc after reconst.	400 µg/ml
Additional dilutions	1:10		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	80 µg/ml	91 µg/ml	91 µg/ml
Cut-off	≥200 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

Order code	Name	C <sub>7</sub> H <sub>5</sub> NaO <sub>3</sub>	
<b>BAG2-C114</b>	<b>Salicylic acid</b>	MW 160.1	
		sodium salt	
Trade Names	Metabolite of Aspirin. Alysine; Idocyl novum; Enterosalicyl; Magsalyl; Parbocyl-Rev; Sodium 2-hydroxybenzoate; 2-Hydroxybenzoic acid		
Content per vial	250 µg	Conc after reconst.	1 mg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	200 µg/ml	227 µg/ml	227 µg/ml
Cut-off	≥ 120 pg/ml	≥5% CD63; SI*≥2	≥10% CD63+CD203c

Order code	Name		
<b>BAG2-C1</b>	<b>Penicillin G</b>		$C_{16}H_{17}N_2NaO_4S$ MW 356.4 sodium salt

Trade Names	Benzylpenicillin Sodium; Sodium Penicillin II; Sodium Benzylpenicillin; Pencillinic Acid Sodium Salt; Penicillin: American Penicillin; Nalpen G; Novocillin; Penilaryn; Pen-A-Brasive; Veticillin		
Content per vial	1 mg	Conc after reconst.	4 mg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	800 µg/ml	910 µg/ml	910 µg/ml
Cut-off	≥ 50 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

Order code	Name		
<b>BAG2-C11</b>	<b>Benzylpenicilloyl-Polylysine (PPL)</b>		

Trade Names	Penicilloypolylysine; PPL; Cilligen; Pre-Pen; Testarpen. The BPO (Benzylpenicilloyl)-group is the major determinant of Typ I allergies against Benzylpenicillin and related antibiotics		
Content per vial	12.5 µg	Conc after reconst.	50 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	10 µg/ml	11.5 µg/ml	11.5 µg/ml
Cut-off	≥50 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

Order code	Name		
<b>BAG2-C12</b>	<b>Minor Determinant Mix (MDM)</b>		

Trade Names	MDM; Mixture of Sodium-Benzylpenicillin and di-Sodium-Benzylpenicilloat		
Content per vial	250 µg	Conc after reconst.	1 mg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	200 µg/ml	225 µg/ml	225 µg/ml
Cut-off	≥ 40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{16}H_{18}N_2O_5SK$
<b>BAG2-C2</b>	<b>Penicillin V</b>	MW 388.5 potassium salt

Trade Names	Acipen-V; Apopen; Distaquaine V; Eskacillin V; Fenacillin; Fenospin; Meropenin; Orocillin; Oratren; Oспен; Pen-Oral; Pen-Vee; Phenopenicillin; P-Mega-Tablinen; Stabicillin; V-Gil; V-Gillin; Vebecillin		
Content per vial	1 mg	Conc after reconst.	4 mg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	800 µg/ml	910 µg/ml	910 µg/ml
Cut-off	≥ 40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{16}H_{19}N_3O_8SZn$
<b>BAG2-C3</b>	<b>Cephalosporin C</b>	MW 478.8 zinc salt

Trade Names	None. Cephalosporin C is the basic structure of many antibiotics such as Cefazolin and Cefamandole.		
Content per vial	25 µg	Conc after reconst.	100 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	20 µg/ml	22.5 µg/ml	22.5 µg/ml
Cut-off	≥ 40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{19}H_{17}N_6O_5S_2$
<b>BAG2-C31</b>	<b>Cefamandole</b>	MW 462.5 formate sodium salt

Trade Names	Bergacef; Cedol; Cefam; Cefiran; Cemando; Cemandil; Fado; Kefadol; Kefandol; Lampomandol; Mandokef; Mandol; Mandolsan; Neocefal; Pavecef		
Content per vial	625 µg	Conc after reconst.	2.5 mg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	500 µg/ml	570 µg/ml	570 µg/ml
Cut-off	≥ 80 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

Order code	Name			$C_{14}H_{14}N_6O_4S_3$
<b>BAG2-C32</b>	<b>Cefazoline</b>			MW 454.5
				sodium salt
Trade Names	Acef; Ancef; Atirin; Biazolina; Bor-Cefazol; Cefacidal; Cefamedin; Cefamezin; Cefazil; Cefazina; Elzogram; Firmacef; Gramaxin; Kefzol; Lampocef; Liviclina; Totacef; Zolicef			
Content per vial	625 µg	Conc after reconst.	2.5 mg/ml	
Additional dilutions	1:5			
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>	
Conc. in stimulation	500 µg/ml	570 µg/ml	570 µg/ml	
Cut-off	≥ 80 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2	

Order code	Name			$C_{16}H_{16}N_4O_8S$
<b>BAG2-C33</b>	<b>Cefuroxime</b>			MW 424.4
				sodium salt
Trade Names	Anaptivan; Biociclin; Biofurex; Bioxima; Cefamar; Cefoprom; Ceftin; Cefumax; Cefurax; Cefurex; Cefurin; Cepacine; Curocef; Curoxim; Duxima; Elobact; Kefurox; Novocef; Spectazole; Zinacef; Zinat; Zinnat			
Content per vial	625 µg	Conc after reconst.	2.5 mg/ml	
Additional dilutions	1:5			
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>	
Conc. in stimulation	500 µg/ml	570 µg/ml	570 µg/ml	
Cut-off	≥ 40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2	

Order code	Name			$C_{15}H_{14}ClN_3O_4S$
<b>BAG2-C34</b>	<b>Cefaclor</b>			MW 367.8
Trade Names	Alenfral; Alfacet; Alfatil; Ceclor; Distaclor; Kefral; Keflor; Panacef; Panoral; Raniclор			
Content per vial	0.75 mg	Conc after reconst.	3 mg/ml	
Additional dilutions	1:10			
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>	
Conc. in stimulation	600 µg/ml	682 µg/ml	682 µg/ml	
Cut-off	≥ 40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2	

<b>Order code</b>	<b>Name</b>	$C_{18}H_{16}N_8Na_2O_7S_3$
<b>BAG2-C35</b>	<b>Ceftriaxone</b>	MW 598.6
		sodium salt

Trade Names Amplospec; Antibacin; Cefotrix; Rocefin; Rocephin; Tercefon

Content per vial	1 mg	Conc after reconst.	4 mg/ml
Additional dilutions	1:10		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	720 µg/ml	818 µg/ml	818 µg/ml
Cut-off	≥ 40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{18}H_{20}FN_3O_4$
<b>BAG2-C36</b>	<b>Levofloxacin</b>	MW 361.4

Trade Names Cravit, Elequine, Floxel, Iquix, Leroxacin, Lesacin, Levaquin, Levokacin, Levox, Levoxacin, Mosardal, Nofaxin, Oftaquix, Quixin, Reskuin, Tavanic, Volequin

Content per vial	200 µg	Conc after reconst.	800 µg/ml
Additional dilutions	1:10		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	160 µg/ml	182 µg/ml	182 µg/ml
Cut-off	≥40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{10}H_{11}N_3O_3S$
<b>BAG2-C61</b>	<b>Sulfamethoxazole</b>	MW 253.3

Trade Names Gantanol; Sinomin. Mixture with Trimethoprim: Bactrim; Berlocid; Cotrim; Cotrimox; Cotrimstada; Cotrimoxazol; Drylin; Escoprim; Eusaprim; Kepinol; Lagatrim; Nopil; Sigaprim; TMS forte

Content per vial	25 µg	Conc after reconst.	100 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	20 µg/ml	22.5 µg/ml	22.5 µg/ml
Cut-off	≥ 50 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

Order code	Name			$C_{14}H_{18}N_4O_3$ MW 290.3
<b>BAG2-C62</b>	<b>Trimethoprim</b>			
Trade Names	Infectotrimet; Motrim; Triprim. Mixture with Sulfamethoxazole: Bactrim; Berlocid; Cotrim; Cotrimox; Cotrimstada; Cotrimoxazol; Drylin; Escoprim; Eusaprim; Kepinol; Lagatrim; Nopil; Sigaprim; TMS forte			
Content per vial	25 µg	Conc after reconst.	100 µg/ml	
Additional dilutions	1:5			
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>	
Conc. in stimulation	20 µg/ml	22.7 µg/ml	22.7 µg/ml	
Cut-off	≥ 40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2	

Order code	Name			$C_{22}H_{24}N_2O_8$ MW 444.4
<b>BAG2-C75</b>	<b>Tetracycline</b>			
Trade Names	Actisite-Dental; Fluorex Plus; Imex; Mysteclin; Sumycin; Tefilin			
Content per vial	25 µg	Conc after reconst.	100 µg/ml	
Additional dilutions	1:5			
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>	
Conc. in stimulation	20 µg/ml	22.7 µg/ml	22.7 µg/ml	
Cut-off	≥ 90 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2	

Order code	Name			$C_{17}H_{18}FN_3O_3$ MW 331.4
<b>BAG2-C81</b>	<b>Ciprofloxacin</b>			
Trade Names	Baycip; Ciflox; Giloxan; Ciprinol; Cipro; Ciprobay; Ciproxan; Ciproxin; Flociprin; Septicide; Velmonit			
Content per vial	25 µg	Conc after reconst.	100 µg/ml	
Additional dilutions	1:5			
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>	
Conc. in stimulation	20 µg/ml	22.5 µg/ml	22.5 µg/ml	
Cut-off	≥ 90 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2	

<b>Order code</b>	<b>Name</b>			$C_{16}H_{19}N_3O_4S$
<b>BAG2-C203</b>	<b>Ampicillin</b>			MW 388.5

Trade Names	Albipen; Amblosin; Ampicin; Ampilar; Ampitab; Amplital; Austrapen; Britacil; Cetampin; Doktacillin; Omnipen; Penbristol; Penbritin; Pentrex; Pentrexyl; Princillin; Principen; Tokiocillin; Totacillin; Ukapen; Vidopen			
Content per vial	2.5 mg	Conc after reconst.	10 mg/ml	
Additional dilutions	1:5			
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>	
Conc. in stimulation	2 mg/ml	2.27 mg/ml	2.27 mg/ml	
Cut-off	≥ 70 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2	

<b>Order code</b>	<b>Name</b>			$C_{16}H_{19}N_3O_5S$
<b>BAG2-C204</b>	<b>Amoxicillin</b>			MW 365.4

Trade Names	Alfamox; Amocilline; Amoxi; Amoxidal; Ardine; Betamox; Clamoxyl; Dura X; Flemoxin; Hiconcil; Ibiamax; Moxal; Neamoxyl; Ospamox; Pasetocin; Penimox; Raylina; Sawacillin; Silamox; Trimox; Ultimox; Zamocillin			
Content per vial	625 µg	Conc after reconst.	2.5 mg/ml	
Additional dilutions	1:5			
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>	
Conc. in stimulation	500 µg/ml	570 µg/ml	570 µg/ml	
Cut-off	≥ 100 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2	

<b>Order code</b>	<b>Name</b>			$C_{38}H_{69}NO_{13}$
<b>BAG2-CCLA</b>	<b>Clarithromycin</b>			MW 748.0

Trade Names	Biaxin; Clathromycin; Cyllind; Klacid; Klaricid; Macladin; Mavid; Naxy; Veclam; Zeclar			
Content per vial	25 µg	Conc after reconst.	100 µg/ml	
Additional dilutions	1:5			
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>	
Conc. in stimulation	20 µg/ml	22.7 µg/ml	22.7 µg/ml	
Cut-off	≥70 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2	

# Antibiotics

**Order code**      **Name**  
**BAG2-CCLAV**      **Clavulanic acid**

$C_8H_9NO_5$   
MW 199.16  
potassium salt

Trade Names                  Combined with amoxicillin (co-amoxiclav, trade names Augmentin, Synulox, and others) or ticarcillin

Content per vial            0.625 mg                                  Conc after reconst.                  2.5 mg/ml

Additional dilutions      1:5

Assay                          **CAST® ELISA**                          **Flow CAST®**                          **Flow CAST® highsens**

Conc. in stimulation      500 µg/ml                                  570 µg/ml                                  570 µg/ml

Cut-off                          ≥40 pg/ml                                  ≥15% CD63                                  ≥15% CD63+CD203c

**Order code**      **Name**  
**BAG2-CRIF**      **Rifampicin**

$C_{43}H_{58}N_4O_{12}$   
MW 823.0

Trade Names                  Rifampicin; Rifaldazine, Rifamycin AMP; R/AMP; Abrifam; Eremfat; Rifa; Rifadine; Rifaldine Rifaldin, Rifaprodin; Rifoldin; Rimactane

Content per vial            25 µg    Conc after reconst.                  100 µg/ml

Additional dilutions      1:5

Assay                          **CAST® ELISA**                          **Flow CAST®**                          **Flow CAST® highsens**

Conc. in stimulation      20 µg/ml    22.7 µg/ml    22.7 µg/ml

Cut-off                          ≥ 100 pg/ml                                  ≥5% CD63; SI\*≥2                                  ≥5% CD63+CD203c; SI\*≥2

# Antiseptics

**Order code**      **Name**

**BAG2-CCHX**      **Chlorhexidine**

$C_{22}H_{30}Cl_2N_{10}x_2(C_6H_{12}O_7)$   
MW 897.8

Digluconate

**Trade Names**      Atoseptal; Chlorhexamed; Corsodyl; Dentohexin; Dosisseptine; Hexamedal; Hibidil; Hibiscrub; Lifo-Scrub; Periochip

<b>Content per vial</b>	1.25 µg	<b>Conc after reconst.</b>	5 µg/ml
-------------------------	---------	----------------------------	---------

**Additional dilutions**    1:5

<b>Assay</b>	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
--------------	--------------------	-------------------	----------------------------

<b>Conc. in stimulation</b>	1 µg/ml	1.14 µg/ml	1.14 µg/ml
-----------------------------	---------	------------	------------

<b>Cut-off</b>	≥100 pg/ml	≥5% CD63; SI*≥2	≥10% CD63+CD203c
----------------	------------	-----------------	------------------

# Contrast Media

<b>Order code</b>	<b>Name</b>			$C_{20}H_{28}I_3N_3O_9$
<b>BAG2-CIOBI</b>	<b>Iobitridol</b>			MW 835.2

Trade Names            Xenetix

Content per vial	1.25 mg	Conc after reconst.	5 mg/ml
------------------	---------	---------------------	---------

Additional dilutions

Assay	CAST® ELISA	Flow CAST®	Flow CAST® highsens
Conc. in stimulation	1 mg/ml	1.1 mg/ml	1.1 mg/ml
Cut-off	≥50 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>			$C_{35}H_{44}I_6N_6O_{15}$
<b>BAG2-CIODI</b>	<b>Iodixanol</b>			MW 1550.2

Trade Names            Visipaque

Content per vial	1.25 mg	Conc after reconst.	5 mg/ml
------------------	---------	---------------------	---------

Additional dilutions

Assay	CAST® ELISA	Flow CAST®	Flow CAST® highsens
Conc. in stimulation	1 mg/ml	1.1 mg/ml	1.1 mg/ml
Cut-off	≥200 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>			$C_{19}H_{26}I_3N_3O_9$
<b>BAG2-CIOHE</b>	<b>Iohexol</b>			MW 821.1

Trade Names            Omnipaque; Accupaque

Content per vial	1.25 mg	Conc after reconst.	5 mg/ml
------------------	---------	---------------------	---------

Additional dilutions

Assay	CAST® ELISA	Flow CAST®	Flow CAST® highsens
Conc. in stimulation	1 mg/ml	1.1 mg/ml	1.1 mg/ml
Cut-off	≥50 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

# Contrast Media

<b>Order code</b>	<b>Name</b>			$C_{17}H_{22}I_3N_3O_8$
<b>BAG2-CIOME</b>	<b>Iomeprol</b>			MW 777.1

Trade Names Iomeron

Content per vial 1.25 mg Conc after reconst. 5 mg/ml

Additional dilutions

Assay	CAST® ELISA	Flow CAST®	Flow CAST® highsens
Conc. in stimulation	1 mg/ml	1.1 mg/ml	1.1 mg/ml
Cut-off	≥50 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>			$C_{17}H_{22}I_3N_3O_8$
<b>BAG2-CIOPA</b>	<b>Iopamidol</b>			MW 777.1

Trade Names Iopamiro; Isovue; Iopamiron; Niopam; Scanlux

Content per vial 1.25 mg Conc after reconst. 5 mg/ml

Additional dilutions

Assay	CAST® ELISA	Flow CAST®	Flow CAST® highsens
Conc. in stimulation	1 mg/ml	1.1 mg/ml	1.1 mg/ml
Cut-off	≥160 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>			$C_{18}H_{24}I_3N_3O_8$
<b>BAG2-CIOPR</b>	<b>Iopromide</b>			MW 791.1

Trade Names Ultravist

Content per vial 1.25 mg Conc after reconst. 5 mg/ml

Additional dilutions

Assay	CAST® ELISA	Flow CAST®	Flow CAST® highsens
Conc. in stimulation	1 mg/ml	1.1 mg/ml	1.1 mg/ml
Cut-off	≥200 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{24}H_{21}I_6N_5O_8$
<b>BAG2-CIOXA</b>	<b>Ioxaglate</b>	MW 1268.9

Trade Names Hexabrix; Ioxaglic acid

Content per vial	1.25 mg	Conc after reconst.	5 mg/ml
------------------	---------	---------------------	---------

Additional dilutions

Assay	CAST® ELISA	Flow CAST®	Flow CAST® highsens
Conc. in stimulation	1 mg/ml	1.1 mg/ml	1.1 mg/ml
Cut-off	≥50 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{12}H_{11}I_3N_2O_5$
<b>BAG2-CIOXI</b>	<b>Ioxitalamate</b>	MW 643.9

Trade Names Telebrix; Ioxitalamic acid

Content per vial	1.25 mg	Conc after reconst.	5 mg/ml
------------------	---------	---------------------	---------

Additional dilutions

Assay	CAST® ELISA	Flow CAST®	Flow CAST® highsens
Conc. in stimulation	1 mg/ml	1.1 mg/ml	1.1 mg/ml
Cut-off	≥200 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

# Local Anesthetics

Order code	Name			$C_{13}H_{20}N_2O_3 \times HCl$
<b>BAG2-CART</b>	<b>Articaine</b>			MW 320.4
hydrochloride				
Trade Names	Astracaine; Rudocain; Septanest; Septocaine; Ubistesin; Ultracaine			
Content per vial	100 µg	Conc after reconst.	400 µg/ml	
Additional dilutions	1:5			
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>	
Conc. in stimulation	80 µg/ml	91 µg/ml	91 µg/ml	
Cut-off	≥ 40 pg/ml	≥5% CD63; SI*≥2	≥10% CD63+CD203c	

Order code	Name			$C_{18}H_{28}N_2O \times HCl$
<b>BAG2-CBUP</b>	<b>Bupivacaine</b>			MW 324.9
hydrochloride				
Trade Names	Carbostesin; Marcain; Marcaine; Sensorcaine; Vivacaine.			
Content per vial	100 µg	Conc after reconst.	400 µg/ml	
Additional dilutions	1:10			
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>	
Conc. in stimulation	80 µg/ml	91 µg/ml	91 µg/ml	
Cut-off	≥ 40 pg/ml	≥5% CD63; SI*≥2	≥10% CD63+CD203c	

Order code	Name			$C_{15}H_{22}N_2O \times HCl$
<b>BAG2-CMEP</b>	<b>Mepivacaine</b>			MW 282.8
hydrochloride				
Trade Names	Carbocaina; Isogaine; Meaverin; Mecain; Mepibil; Mepicaton; Mepident; Mepidont; Mepiforan; Mepihexal; Mepimynol; Mepivacaina; Mepivamol; Mepivastesin; Scandicain; Scandicaine; Tevacaine			
Content per vial	100 µg	Conc after reconst.	400 µg/ml	
Additional dilutions	1:10			
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>	
Conc. in stimulation	80 µg/ml	91 µg/ml	91 µg/ml	
Cut-off	≥40 pg/ml	≥5% CD63; SI*≥2	≥10% CD63+CD203c	

<b>Order code</b>	<b>Name</b>
<b>BAG2-CLID</b>	<b>Lidocaine</b>

$C_{14}H_{20}N_2O Na_2$
MW 278.3
sodium salt

Trade Names	Cuivasil; Duncaine; Leostesin; Lidesthesin, Lidostesin; Lignavet; Odontalg; Rucaina; Sedagul; Xylocaine; Xylocard; Xylocitin; Xyloneural; Xylotox;		
-------------	--	--	--

Content per vial	62.5 µg	Conc after reconst.	250 µg/ml
------------------	---------	---------------------	-----------

Additional dilutions	1:5
----------------------	-----

Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
-------	--------------------	-------------------	----------------------------

Conc. in stimulation	50 µg/ml	57 µg/ml	57 µg/ml
----------------------	----------	----------	----------

Cut-off	≥40 pg/ml	≥5% CD63; SI*≥2	≥10% CD63+CD203c
---------	-----------	-----------------	------------------

# Myorelaxants

<b>Order code</b>	<b>Name</b>	$C_{65}H_{82}N_2O_{18}S_2$
<b>BAG2-CATR</b>	<b>Atracurium</b>	MW 1243.5
		besylate salt

Trade Names Acuron; Aratan; Atratan; Atrelax; Atrican; Tracrium; Trarium; Wellcome 33-A-74

Content per vial	1.25 µg	Conc after reconst.	5 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	1 µg/ml	1.2 µg/ml	1.2 µg/ml
Cut-off	≥50 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{53}H_{72}N_2O_{12} \times 2 C_6H_5SO_3$
<b>BAG2-CCAT</b>	<b>Cis-Atracurium</b>	MW 929.1
		Besylate

Trade Names Nimbex

Content per vial	100 µg	Conc after reconst.	400 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	80 µg/ml	91 µg/ml	91 µg/ml
Cut-off	≥ 40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{58}H_{80}Cl_2N_2O_{14}$
<b>BAG2-CMIV</b>	<b>Mivacurium</b>	MW 1100.2

Trade Names Mivacron

Content per vial	250 µg	Conc after reconst.	1 mg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	200 µg/ml	227 µg/ml	227 µg/ml
Cut-off	≥40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

# Myorelaxants

<b>Order code</b>	<b>Name</b>			$C_{35}H_{60}Br_2N_2O_4$
<b>BAG2-CPAN</b>	<b>Pancuronium</b>			MW 732.7

Trade Names            Mioblock; Pancron; Pancur; Pancuronium; Parium; Pavulon

Content per vial	250 µg	Conc after reconst.	1 mg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	200 µg/ml	227 µg/ml	227 µg/ml
Cut-off	≥ 110 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>			$C_{12}H_{18}O$
<b>BAG2-CPRO</b>	<b>Propofol</b>			MW 178.3

Trade Names            Diprivan; Disoprivan; Fresofol; Pofol; Propofol; Rapinovet; Recofol

Content per vial	250 µg	Conc after reconst.	1 mg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	200 µg/ml	227 µg/ml	227 µg/ml
Cut-off	≥50 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>			$C_{32}H_{53}N_2O_4$
<b>BAG2-CROC</b>	<b>Rocuronium</b>			MW 371.4

Trade Names            Esmeron; Zemuron

Content per vial	250 µg	Conc after reconst.	1 mg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	200 µg/ml	227 µg/ml	227 µg/ml
Cut-off	≥70 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>
<b>BAG2-CSUX</b>	<b>Suxamethonium</b>

$C_{14}H_{30}Cl_2N_2O_4$   
MW361

chloride salt

Trade Names            Lysthenon; Midarine; Pantolax; Succinolin

Content per vial	2.50 mg	Conc after reconst.	10 mg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	2 mg/ml	2.27 mg/ml	2.27 mg/ml
Cut-off	≥40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>
<b>BAG2-CVEC</b>	<b>Vecuronium</b>

$C_{34}H_{57}BrN_2O_4$   
MW 637.7

bromid salt

Trade Names            Musculax; Norcuron

Content per vial	62.5 µg	Conc after reconst.	250 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	50 µg/ml	51 µg/ml	51 µg/ml
Cut-off	≥40 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>
<b>BAG2-CNBX</b>	<b>Neuromuscular Blockers Mix</b>

Ingredients	Content per vial µg	Con. after reconst. µg/ml
Suxamethonium	795	3180
Recuronium	43	171
Vecuronium	96	383
Pancuronium	95	381
Atracurium	12	50
Mivacurium	10	41
Cisatracurium	20	82

Content per vial	1.072 mg	Conc after reconst.	4.3 mg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	0.86 mg/ml	0.97 mg/ml	0.97 mg/ml
Cut-off	≥ 200 pg/ml	≥5% CD63; SI*≥2	≥10% CD63+CD203c

# Proton Pump Inhibitors

<b>Order code</b>	<b>Name</b>	$C_{17}H_{19}N_3O_3S$
<b>BAG2-COME</b>	<b>Omeprazole</b>	MW 345.4

Trade Names	Antra, Gastroloc, Losec, Losectil, Lozeprel, Mopral, Omepral, Omez, Prilosec, Ulcozol		
Content per vial	10 µg	Conc after reconst.	40 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	8 µg/ml	9.1 µg/ml	9.1 µg/ml
Cut-off	140 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{34}H_{36}MgN_6O_6S_2$
<b>BAG2-CESO</b>	<b>Esomeprazole</b>	MW 713.1
		magnesium hydrate

Trade Names	Nexium, Nexiam, Nexium IV, Axagon, Esopral, Lucen		
Content per vial	10 µg	Conc after reconst.	40 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	8 µg/ml	9.1 µg/ml	9.1 µg/ml
Cut-off	100 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

<b>Order code</b>	<b>Name</b>	$C_{16}H_{14}F_2N_3O_4S \times Na$
<b>BAG2-CPANT</b>	<b>Pantoprazole</b>	MW 405.4
		sodium salt

Trade Names	Astropan, Controloc, Inipomp, Pantomed, Pantecta, Pantoloc, Pantozol, Pantec, Pantoprazole, Pantor, Protium, Protonix, Protonix IV, Somac, Tectac,		
Content per vial	50 µg	Conc after reconst.	200 µg/ml
Additional dilutions	1:5		
Assay	<b>CAST® ELISA</b>	<b>Flow CAST®</b>	<b>Flow CAST® highsens</b>
Conc. in stimulation	40 µg/ml	45.5 µg/ml	45.5 µg/ml
Cut-off	110 pg/ml	≥5% CD63; SI*≥2	≥5% CD63+CD203c; SI*≥2

# Proton Pump Inhibitors

**Order code**      **Name**  
**BAG2-CLAN**      **Lansoprazole**

$C_{16}H_{14}F_2N_3O_2S$   
 MW 396.4

Trade Names      Aگوpton, Lanciprol, Lansobene, Lansoloc, Lansoprazol, Lansox, Lanton, Lanzo, Lanzopral, Lanzor, Lanzostad, Limpidex, Monolimum, Prevacid, Zoton, Zomel

Content per vial      50 µg      Conc after reconst.      200 µg/ml

Additional dilutions      1:5

Assay      **CAST® ELISA**      **Flow CAST®**      **Flow CAST® highsens**

Conc. in stimulation      40 µg/ml      45.5 µg/ml      45.5 µg/ml

Cut-off      200 pg/ml      ≥5% CD63; SI\*≥2      ≥5% CD63+CD203c; SI\*≥2

**Order code**      **Name**  
**BAG2-CRAB**      **Rabeprazole**

$C_{18}H_{20}N_3O_3S \times Na$   
 MW 381.4  
 sodium salt

Trade Names      AcipHex, Pariet, Rebezole, Rabezole

Content per vial      50 µg      Conc after reconst.      200 µg/ml

Additional dilutions      1:5

Assay      **CAST® ELISA**      **Flow CAST®**      **Flow CAST® highsens**

Conc. in stimulation      40 µg/ml      45.5 µg/ml      45.5 µg/ml

Cut-off      100 pg/ml      ≥5% CD63; SI\*≥2      ≥5% CD63+CD203c; SI\*≥2

BÜHLMANN Laboratories AG  
Baselstrasse 55  
CH-4124 Schönenbuch/Basel  
Switzerland

Phone +41 61 487 12 12  
Fax orders +41 61 487 12 99  
info@buhlmannlabs.ch  
www.buhlmannlabs.ch

**Ordering Code:**

Flow CAST® FK-CCR 100 tests  
Flow CAST® highsens FK-HSBAT 100 tests

CAST® ELISA EK-CAST 96 wells  
EK-CAST5 480 wells

CAST® is a registered trademark of  
BÜHLMANN Laboratories AG