



Distributor Presentation


December 2023

Updated by George Skinner

Internal, Privileged and Confidential
Not for distribution

PRODUCT

FASTinov workflow: AST from blood cultures in 2 hours



FASTinov AST kits are CE-marked, approved for sale in the European Union under IVD Directive. They are not FDA-cleared and are not yet available for sale in the U.S. Availability in each country depends on local regulatory clearance.


Kit configuration – current design



Documentation



Secondary box



Tube A with Hemolytic agent X100



1x vials with flat base
w/close cap 0.5ml
X100


In paper box, labelled



AST panels*




8x FAST panels (v1, commercial)
closed and with primary label
attached



Plastic cuvettes



8x cuvettes




FASTinov technology: flow cytometry and phenotypic cell lesion study


FASTinov AST panels



Flow cytometry



bioFAST software



FASTinov AST kits are CE-marked, approved for sale in the European Union under IVD Directive. They are not FDA-cleared and are not yet available for sale in the U.S. Availability in each country depends on local regulatory clearance.

Antibiotic panel - FASTgramneg

Antibiotics FASTgramneg	<i>Enterobacterales</i>		<i>Pseudomonas spp</i>		<i>Acinetobacter spp</i>	
	EUCAST	CLSI	EUCAST	CLSI	EUCAST	CLSI
Ampicillin	✓	✓	-	-	-	-
Amoxicillin-clavulanic acid	✓	✓	-	-	-	-
Cefotaxime	✓	✓	-	-	-	-
Ceftazidime	✓	✓	✓	✓	-	-
Cefepime	✓	✓	✓	✓	-	-
Piperacillin-tazobactam	✓	✓	✓	✓	-	✓
Ceftalozane-tazobactam	✓	✓	✓	✓	-	-
Ceftazidime-avibactam	✓	✓	✓	✓	-	-
Meropenem	✓	✓	-	-	-	-
Ciprofloxacin	✓	✓	✓	✓	✓	✓
Gentamicin	✓	✓	-	✓	✓	✓
Amikacin	✓	✓	✓	✓	✓	✓

- Provides AST for **12 antibiotics** for *Enterobacterales*, *Pseudomonas* and *Acinetobacter* spp.
- Detects main mechanism of resistance ESBL for EB Group 1
- Screens for presence of ESBL (for EB Group 2), pAmpC and carbapenemases



Antibiotic panel - FASTgrampos

Antibiotics FASTgrampos	Staphylococcus spp		Enterococcus spp	
	EUCAST	CLSI	EUCAST	CLSI
Penicillin*	✓	✓	✓	✓
Ampicillin	-	-	✓	✓
Cefoxitin**	✓	✓	-	-
Oxacillin***	✓	✓	-	-
Vancomycin****#	✓	✓	✓	✓
Linezolid	✓	✓	✓	✓
Gentamicin	✓	✓	-	-
Gentamicin high level#	-	-	✓	✓

* Only for *S. aureus* according to EUCAST;

** Not applicable to *S. epidermidis*;

*** Only for *S. epidermidis*; **** MIC value for *S. aureus*;

Only for *E. faecalis*

- Provides AST for **7 antibiotics** for *Staphylococcus* and *Enterococcus* spp.
- Provides **MIC** for **Vancomycin** for *S. aureus*

Test accuracy



	EUCAST							CLSI						
	N	EA%	CA%	mE	ME	VME	Reproducibility	N	EA%	CA%	mE	ME	VME	Reproducibility
FASTgramneg	2570	NA	98.9	0.0%	1.2%	0.9%	99.5%	2677	NA	97.9	1.3%	0.8%	0.7%	99.3%
FASTgrampos	857	91.5	97.8	0.2%	2.8%	0.4%	99.4%	931	91.5	97.9	0.3%	2.6%	0.3%	99.8%

n- number of strains

EA (%) - Percentage of Essential agreement (only for colistin and vancomycin)

CA (%) - Percentage of Categorical Agreement

mE- minor Error; ME- Major Error; VME- Very Major Error

NA- not applicable


		Detection of mechanisms of resistance			
		N	Sensitivity	Specificity	Accuracy
Screening	ESBL(EB group I)	45	95.7%	100%	99.3%
	ESBL(EB group II)	17	100%	100%	100%
	pAmpC	37	100%	100%	100%
	Carbapenemases	52	92.2%	95.1%	94.1%

Source: IFU


FASTinov AST kits are CE-marked, approved for sale in the European Union under IVD Directive. They are not FDA-cleared and are not yet available for sale in the U.S. Availability in each country depends on local regulatory clearance.

Space needed

Example A




Example B



POSITIONING & PRICING


Competitive positioning – rapid AST



Competitor messages – Jul 2023



PHENOTYPIC AST IN 2H
FROM BLOOD CULTURES



Every
competitor
wants to claim
"rapid result"
ability



Designed to save lifetimes

Fast and comprehensive – when time is life

- Fully automated
- True MIC results
- Extensive AST capabilities



VITEK® REVEAL™

RAPID AST FOR BLOODSTREAM INFECTIONS

VITEK® REVEAL™ integrated within the bioMérieux Sepsis Solution enables same-day Antimicrobial Susceptibility Testing (AST) results to help you manage bloodstream infections:

- Delivers phenotypic AST results in an average of 5.5 hours¹
- Directly from positive blood culture
- Wide antimicrobial coverage for Gram negative bloodstream infections



It's okay to speed.
You're providing
MIC results.

Fast

Fully Automated

Actionable results



dRAST™, direct & Rapid Antimicrobial Susceptibility Test

dRAST provides MIC-based phenotypic antimicrobial susceptibility testing direct from positive blood culture in as low as 4 hours; reducing time to results by up to 2 days compared to conventional method. The rapid results allow for earlier optimization of antibiotics in critically-ill patients with bloodstream infections and sepsis.

- ✓ Dedicated to positive blood culture samples
- ✓ Provides phenotypic MIC in as low as 4 hours
- ✓ Random access with up to 12 samples simultaneously (15 panel positions)
- ✓ Expert system on board with choice of guidelines
- ✓ Easy start with no McFarland required, no sample prep
- ✓ 2 panels: 1 Gram Neg. + 1 Gram Pos.
- ✓ Easy to use, fast to operate
- ✓ No daily maintenance

Speed

Speed

MIC

Speed

MIC

Automated

Speed

Speed

MIC

Automated

Speed

MIC

Automated

Competitive positioning: rapid AST competitors



	FASTinov	Gradientech	Qlinea	SpecificDx	Accelerate Dx	Quantamatrix
Time to AST results	2h SUPERIORITY	2-4h ADVANTAGE	6-8h WEAKNESS	6-8 h WEAKNESS	6-8h WEAKNESS	6-8h WEAKNESS
Coverage of antibiotics panel	GN: 12 Abx GP: 7 Abx STRENGTH	GN: 13 Abx GP: - WEAKNESS	GN: 23 Abx GP: -- WEAKNESS	GN: Available GP: WEAKNESS	GN: 12 Abx GP: 5 Abx STRENGTH	GN: 12 Abx GP: 18 Abx STRENGTH
Throughput (1 instrument)	1 at a time 10-12 AST in 8h STRENGTH	1 at a time 10-12 AST in 8h STRENGTH	4 simultaneous 4 AST in 8h STRENGTH	4 simultaneous 4 AST in 8h WEAKNESS	1 at a time 1 AST in 8h WEAKNESS	16 simultaneous 16 AST in 8h STRENGTH
Usability	5-10' hands-on time PARITY	2'30" hands-on time STRENGTH	2' hands-on time STRENGTH	2'30" hands-on time STRENGTH	2' hands-on time STRENGTH	2' hands-on time STRENGTH

New category:
"Ultra Rapid"


"rapid AST" category

AST with TTR 2h =>
Same day


"Rapid" AST with TTR 6-8 hours =>
Most results delivered **next day**

Competitive positioning - accuracy


Gram-negative



Gram-positive



Total



Comparison of antibiotic panel – gram-negative



			FASTinov	Gradientech	QLinea	Accelerate	Reveal	Quantamatrix
Penicillins		Ampicillin	✓		✓		✓	✓
		Piperacillin					✓	✓
Cephalosporins	1st gen cephalosporins	Cefazolin			✓		✓	
	2nd gen cephalosporins	Cefuroxime			✓		✓	
		Cefoxitin					✓	
	3rd gen cephalosporins	Cefotaxime	✓	✓	✓		✓	✓
		Ceftazidime	✓	✓	✓	✓	✓	✓
		Ceftriaxone			✓	✓		
	4th gen cephalosporins	Cefepime	✓	✓	✓	✓	✓	✓
Penicillins+Beta Lactamase Inhibitor		Piperacillin/tazobactam	✓	✓	✓	✓	✓	✓
		Amoxicillin/clavulanic acid	✓		✓		✓	✓
		Ampicillin-Sulbactam				✓	✓	✓
Cephalosporin+Beta Lactamase Inhibitor		Ceftolozane/tazobactam	✓		✓		✓	
		Ceftazidime/avibactam	✓	✓	✓		✓	✓
		Ceftazidime/Clavulanic acid	✓				✓	✓
Carbapenems		Meropenem	✓	✓	✓		✓	✓
		Ertapenem			✓	✓	✓	
		Imipenem					✓	✓
Quinolones		Ciprofloxacin	✓	✓	✓	✓	✓	✓
		Levofloxacin			✓		✓	✓
Aminoglycosides		Gentamicin	✓	✓			✓	✓
		Amikacin	✓	✓		✓	✓	✓
		Tobramycin		✓		✓	✓	
Other	Antifolates	Trimethoprim/Sulfamethoxazole.					✓	✓
	Monobactams	Aztreonam			✓	✓	✓	
	Tetracycline	Tigecycline		✓			✓	
	Polymyxin	Colistin		✓				✓

2h

2-4h

6-8h

6-8h

6-8h


6-8h

Comparison of antibiotic panel – gram-positive




Gram-positive			FASTinov	Gradientech	QLinea	Accelerate	Reveal	Quantamatrix
Penicillins		Ampicillin	✓			✓		
		Oxacillin	✓					
		Penicillin	✓					
Cephalosporins	2nd gen cephalosporins	Cefoxitin	✓					
	5th gen cephalosporins	Ceftaroline				✓		
Oxazolidinones		Linezolid	✓			✓		
Lipopeptides		Daptomycin				✓		
Glycopeptides		Vancomycin	✓			✓		
Aminoglycosides		Gentamicin	✓					
			2h	2-4h	6-8h	6-8h	6-8h	6-8h

While the Automated AST market grows with intense rivalry, new segment “**AST results that are actionable in the same shift**” is a unique and differentiated, high-value offer





Competitive pricing: FASTinov intends to be price leader

Price per test



CUSTOMER TARGETS / SEGMENTS

Target Customer Segment

- Large and mid-size Hospitals with clinical pathology Dept.
 - Microbiology/Clinical pathology labs  budget  champion
 - Clinical departments  supporter
 - ICU, E&R, Int Medicine/Infecciology , AMR stewardship, Pediatric ICU
 - Oncology / Transplant departments
 - Purchasing / Administration  gatekeeper
- Group purchasers / Regional purchasers in tender regions  gatekeeper



Attitudinal:

- Wants to improve profile of pathology, especially post-covid
- Wants to be a partner to clinical services
- Believes in the importance of rapid AST to antibiotic stewardship

- Hospitals with rapid identification capability (Maldi-TOF / Mass spectrometry).
- Rapid Molecular ID and susceptibility systems in place indicated that customer has already obtained budget for rapid testing
- Infection control commission
- 500+ beds, or specialized hospitals in oncology, transplants, or pediatric & adult intensive units

Sales funnel – stages

Promotional stage > Demo stage > Commercial stage > Post-selling support

Stage 1

- Relationship building
- Promotion
- Business intelligence
- Tender management

Stage 2

- Evaluation planning
- Installation
- Certification
- Evaluation

Stage 3

- Proposal preparation
- Pricing approval
- Negotiation
- Activation

Stage 4

- Training
- Reorders
- Business intelligence

Gate
1

Evaluation agreed

Gate
2

Negotiation/
tender started

Gate
3

Approval to start

MESSAGING & MARKETING MIX

We have analysed the launch in Europe and we have a couple of important lessons learned.

Rapid antibiotic susceptibility testing (rAST) is a **young field in Europe**.

The Corona pandemic has had continued effects during 2022 with stretched hospital budgets and **limited staff availability**.

Customers want to **"try before you buy"**.

How does the system perform in my clinical setting?

Customers are asking for **clinical and health economic evidence (HEOR)**.

Customers want to **compare the different rAST solutions present in the market**.

Voice of **key opinion leaders** and **scientific evidence** is important.

Q-LINEA 

This view of the market is from a competitor and highlights the situation that we all are facing, quite well.

- For the **critical infections** that require **immediate decisions**
- FASTinov AST provides the **world's fastest AST with results in less than 2h**
- For a **broad panel of species and antibiotics**
- AST results that are **Actionable within the Same Shift**



**PHENOTYPIC AST IN 2H
FROM BLOOD CULTURES**

**ULTRA
RAPID**



www.fastinov.com

**ULTRA
RAPID**



**PHENOTYPIC AST IN 2H
FROM BLOOD CULTURES**

www.fastinov.com

**ULTRA
RAPID**



**PHENOTYPIC AST IN 2H
FROM BLOOD CULTURES**

COMMERCIAL TOOLS

Pricing and discounting scheme – Recommended

Preferred discounting scheme based on **FOC**.

		Price per test, EUR, VAT excluded
List price (per test), EUR		190.00
Instrument discount	0.0%	0.00
Introductory offer discount	5.0%	9.50
Invoice price (per test), EUR		180.50

Instrument placement

2) Availability of cytometer

Cytometers already installed	Discount
Not installed	0%
Installed	10%

Free-of-charge system

Annualized panel consumption*	Nr of FOC tests	Final ASP at min
<= 300	0	181
301 - 600	8	176
601 - 900	32	171
901 - 1200	192	149
1201 - 1500	320	143
> 1500	540	133

Introductory offer: for hospitals starting with FASTinov AST

3) Introductory offer

	Discount
New client	8%
Existing customer	0%

Estimated unit price after FOC, EUR** **149**


corresponding to 192 tests offered free of charge plus 910 tests purchased

Training and certification process

FASTinov provides on-site and off-site certification for operators of the technology.

- Product
- Training and certification
- Software update
- Lab workflow

SOP



Training

Certification

Introduction
presentation

Step-by-step
demo in lab

1st case
observation


Certification



Commercial Start


- Jan-Jun 2023
- 6 demos completed in Portugal
 - POR_Porto_IPO Porto
 - POR_Almada_Garcia de Orta
 - POR_Coimbra_CHUC
 - POR_Vila Real_CHTMAD
 - POR_Amadora_Fernando Fonseca
 - POR_Porto_Santo António
- 139 tests performed

Portugal - Customer demos 2023



IPO Porto

26 Jan → 15 Feb, N=15



19 Apr → 2 May,
N=29



CHTMAD, Vila Real



6 Jun → 20 Jun,
N=12



22 Mar → 12 Apr
N=15



17 May → 2 Jun,
N=30e




HFF, Amadora




23 Feb → 15 Mar, N=49



HGO, Almada



HGO, Almada



CHUC, Coimbra


Events

- 2 – 4 March Figuera de Foz Portugal
- 8-11 March AMCLI Rimini Italy
- 27-30 April 2024 ECCMID Barcelona Spain
- 30 May - 1 June Zaragoza Spain
- ?? Sept Sepsis Update – Weimar, Germany



Online


- **Website**





Field training


- Sharepoint
- Materials
- **Audiovisuals** (Video animation, loop presentation, mini-symposia)
 - Demonstration reports
 - Training and **lab support materials**
 - Website and **online awareness** on the rise



Field training – sharepoint (June 2023)


<https://fastinov.sharepoint.com/sites/FASTinov?OR=Teams-HL&CT=1689327212092&clickparams=eyJBCbHBOYW1lIjoiVGhvbXMtRGVza3RvcCIsIkFwcFZlcnNpb24iOiIiNDE1LzIzMDYwNDExMTYxIiwicGFzRmVkdXJhdGVkVXNlciI6ZmFsc2V9>

Communication of Didactic Resources for Training




1st FASTinov Demonstration

FASTinov team will be present in the laboratory and perform the procedure for the selected samples.




2nd Training

FASTinov team will provide initial training on susceptibility through flow cytometry to the local team.




3rd Hand Over


Local team will be present in the laboratory and perform the procedure for the selected samples.



FASTinov - Fastest workflow f...




FASTinov Training




biofast


Video exemplificativo de como retirar um reporte de suscetibilidade no bioFAST.



PubMed




FASTinov Training



sw-dx


Video exemplificativo de como inserir uma experiência no




PubMed

Evaluation of FASTinov Ultrarapid Flow Cytometry Antimicrobial Susceptibility Testing Directly from...
pubmed.ncbi.nlm.nih.gov

The FASTinov flow cytometry kit, an ultrarapid antimicrobial susceptibility test, was directly evaluated on positive blood cultures (BC) at two sites: (i) FASTinov, S.A., in Porto, Portugal, using BC spiked with well-characterized bacteria, and (ii) ...



FASTinov Training




FASTinov Training

Português (Portugal)


Quiz: Metodologia FASTinov

O objetivo deste formulário é consolidar as informações adquiridas pelos técnicos após formação e experiência na realização do procedimento laboratorial FASTinov.


- Qual é a ordem correta para inserir os códigos de barras de forma a guardar a experiência para leitura no citoemtor de fluxo?
 - ☐ Código do doente Código do painel
 - ☐ Código do painel Código do doente
- Como se devem encontrar os recipientes "Sheath"/"Waste" no início da utilização do citoemtro?
 - ☐ Sheath preenchido até ao limite máximo / Waste vazio
 - ☐ Sheath vazio / Waste preenchido até ao limite máximo
- Qual é o procedimento a adoptar quando o controlo de qualidade falha?
 - ☐ Continuar o procedimento seguinte o "Daily Clean" e voltar a executar o procedimento do controlo de qualidade
 - ☐ Repetir novamente o procedimento do controlo de qualidade
 - ☐ Avançar para o procedimento do "Daily ...



Video exemplificativo de como inserir uma experiência no Citoemtro de fluxo para fazer a sua leitura.




Procedimento para ID por MALDI-TOF



Inicialização do Citoemtro e Manutenção Diária

English (United States)



FASTinov
Next Generation
Antimicrobial Resistance Technologies

FASTinov AST - User Experience Questionnaire

This survey was created to collect the user experience from FASTinov AST test users. The goal is to understand the levels of satisfaction / dissatisfaction to keep improving our product, training, and service.

The survey will take approximately 7 minutes to complete.

**Clinical
marketing**

Bloodstream infection (BSI)

Incidence:

**200-220 episodes / 100,000
population^{1,2}
(growing at 14% per year¹)**

Case fatality rates:

**25.6% at 30 days
46.4% at 1 year²**

Empiric Antibiotic Therapy (EAT)

*N = 5,715 patients
with septic shock³*

**Appropriate
80%**

**Inappropriate
20%**

*“antimicrobials with in
vitro activity for the isolated
pathogen(s) or appropriate
for the underlying clinical
syndrome”*

**Survival
52%**

**Survival
10.3%**

¹ Buetti N, Atkinson A, Marschall J, *et al* Incidence of bloodstream infections: a nationwide surveillance of acute care hospitals in Switzerland 2008–2014 *BMJ Open* 2017;7:e013665. doi: 10.1136/bmjopen-2016-013665

² Schechner, Vered *et al*. One-year mortality and years of potential life lost following bloodstream infection among adults: A nation-wide population based study *The Lancet Regional Health – Europe*, Volume 23, 100511

³ Kumar A, Ellis P, Arabi Y, *et al*. ; Cooperative Antimicrobial Therapy of Septic Shock Database Research Group. Initiation of inappropriate antimicrobial therapy results in a fivefold reduction of survival in human septic shock. *Chest* 2009, 136:1237–48.

Bloodstream infection (BSI)

Incidence:

200-220 episodes / 100,000 population^{1,2}
(growing at 14% per year¹)

Case fatality rates:
25.6% at 30 days
46.4% at 1 year²

**Susceptibility
Testing is
necessary to
inform
appropriate
therapy**

Empiric Antibiotic Therapy (EAT)

*N = 5,715 patients
with septic shock³*

**Appropriate
80%**

*“antimicrobials with in
vitro activity for the isolated
pathogen(s) or appropriate
for the underlying clinical
syndrome”*

**Survival
52%**

**Inappropriate
20%**

**Survival
10.3%**

¹ Buetti N, Atkinson A, Marschall J, *et al* Incidence of bloodstream infections: a nationwide surveillance of acute care hospitals in Switzerland 2008–2014 *BMJ Open* 2017;7:e013665. doi: 10.1136/bmjopen-2016-013665

² Schechner, Vered *et al*. One-year mortality and years of potential life lost following bloodstream infection among adults: A nation-wide population based study *The Lancet Regional Health – Europe*, Volume 23, 100511

³ Kumar A, Ellis P, Arabi Y, *et al*. ; Cooperative Antimicrobial Therapy of Septic Shock Database Research Group. Initiation of inappropriate antimicrobial therapy results in a fivefold reduction of survival in human septic shock. *Chest* 2009, 136:1237–48.

Ultra-rapid susceptibility → involvement in clinical decisions

**Empiric
Antibiotic
Therapy**

Espécie	Klebsiella oxytoca
Painel	FASTgramneg
Tipo de Amostra	Hemocultura
ID Amostra	EXP_20230313_123807137
Número do Lote	22C1122
Técnico	rita

Versão AST	2.1
Incubação Início	13/03/2023 12:10:49
Instrumento Início	13/03/2023 13:17:43
Relatório Final	13/03/2023 14:48:00

Suscetibilidade Antimicrobiana

Fármaco

Ampicilina	R
Amoxicilina/Ácido clavulânico	S
Cefotaxima	S
Ceftazidima	S
Cefepime	S
Piperacilina/Tazobactam	S
Ceftalozane+Tazobactam	S
Ceftazidime+Avibactam	S
Meropenem	S
Ciprofloxacina	S
Gentamicina	S
Amicacina	S

**Optimal
Antibiotic
Therapy**



Thank You!