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


## Literatura

1. Thayer, J.D., and J.E. Martin, Jr. 1966. Improved medium selective for cultivation of *N.gonorrhoeae* and *N. meningitidis*. Public Health Rep. 81:559-562.
2. Martin, J.E., J.H. Armstrong, and P.B. Smith. 1974. New system for cultivation of *Neisseria gonorrhoeae*. Appl. Microbiol. 27:802-805.
3. Martin, J.E., Jr., and J.S. Lewis. 1977. Anisomycin: improved antimycotic activity in modified Thayer- Martin medium. Public Health Lab. 35:53-62.
4. Cross, R.C., M.B. Hoger, R. Neibaur, B. Pasternack, and F.J. Brady. 1971. VCN-inhibited strains of *Neisseria gonorrhoeae*. HSMHA Health Rep. 86:990-992.
5. Phillips, I., D. Humphrey, A. Middleton, and C.S. Nicol. 1972. Diagnosis of gonorrhea by culture on a selective medium containing vancomycin, colistin, nystatin, and trimethoprim (VCNT). A comparison with gram-staining and immunofluorescence. Brit. J. Vener. Dis. 48:287-292.
6. Reichart, C.A., L.M. Rupkey, W.E. Brady, and E.W. Hook III. 1989. Comparison of GC-Lect and modified Thayer-Martin media for isolation of *Neisseria gonorrhoeae*. J. Clin. Microbiol. 27:808-811.
7. Evans, G.L., D.L. Kopyta, and K. Crouse. 1989. New selective medium for the isolation of *Neisseria gonorrhoeae*. J. Clin. Microbiol. 27:2471-2474.
8. Thomson, R.B., and J.M. Miller. 2003. Specimen collection, transport, and processing: bacteriology. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
9. Isenberg, H.D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
10. Janda, W.M., and J.S. Knapp. 2003. *Neisseria* and *Moraxella catarrhalis*. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.

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