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Dátum poslednej revízie: 10.03.2014

**SUMMARY AND EXPLANATION**

Assimilation diagnostic discs are designed for identification of clinically significant species of yeasts.

**PRINCIPLE**

Principle is based on ability of yeasts to assimilate carbohydrates and other chemicals diffused from disc into agar, what causes that grow of strain is stimulated around assimilation disc.

**DIAGNOSTIC KIT CONTAINS**

- tube with diagnostic discs - 500 pcs

Ref.	Name	Description	Ref.	Name	Description
6001	ARA	Utilization of arabinose	6011	RAF	Utilization of raffinose
6002	CEL	Utilization of cellobiose	6012	RHA	Utilization of rhamnose
6003	GAL	Utilization of galactose	6013	RIB	Utilization of ribose
6004	GLU	Utilization of glucose	6014	SOR	Utilization of sorbitol
6005	INO	Utilization of inositol	6015	SUC	Utilization of sucrose
6006	LAC	Utilization of lactose	6016	TRE	Utilization of trehalose
6007	MNS	Utilization of mannose	6017	XYL	Utilization of xylose
6008	MLB	Utilization of melibiose	6018	NIT	Utilization of KNO <sub>3</sub>
6009	MLT	Utilization of maltose	6019	PEP	Utilization of peptone
6010	MLZ	Utilization of melzitose	6020	URE	Utilization of urea

**WARNINGS AND SAFETY PRECAUTIONS**

- For *in vitro* diagnostics use.
- For professional use only.
- Follow the instructions exactly!
- All samples and inoculated products have to be considered as potentially infectious and obvious safety valid regulations must be respected while handling.
- Do not use product out of expiration
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of overexpired or devaluated product liquidate according to valid regulations. Packaging materials are recyclable.

**STORAGE**

Keep diagnostic discs in dark at temperatures between +2 and +25 °C.

**REQUIRED MATERIALS**

- C-agar for auxanogram for discs with ref. 6001-6017, N-agar for discs with ref. 6018-6020
- sterile Petri dish, sterile inoculation loop, water bath
- saline 1 ml

**RECOMMENDED PROCEDURE**

1. Take one full inoculation loop of yeast mass from 24 hours old culture and suspend it in tube with sterile saline (1 ml). Suspensions must be prepared immediately before its usage.
2. Mix prepared suspension with agar for auxanogram, which was melted before, in water bath (t < 52°C). Pour agar on sterile Petri dish, than, after solidification, put discs on agar (max. 7 on Petri dish). Press discs gently, so that whole surface of discs is in contact with medium. Incubate it in thermostat for 24 h at temperatures (+30 ± 2)°C.

**EVALUATION AND INTERPRETATION**

As positive result is interpreted grow of colonies around diagnostic discs seen as round zone. In case of negative result grow of colonies around disc is not visible.

Functionality of discs can be verified on control strains (C.S.):

1. *Candida guilliermondii* CCM 8321 / ATCC 90877
2. *Cryptococcus neoformans* CCM 8312 / ATCC 90112
3. *Klebsiella pneumoniae* CCM 5852 / ATCC 13882

C.S.	A	R	E	G	G	I	L	M	M	M	R	R	R	S	S	T	X	N	P	U		
C.S.	A	R	E	A	L	U	O	C	S	B	T	Z	F	A	B	R	C	E	L	T	P	E
1.	+	+	+	+	-	-	+	+	+	+	+	+	+	v	-	+	+	+	+	-	+	-
2.	-	-	+	+	+	-	+	-	+	+	+	+	v	+	+	+	+	+	+	-	+	+
3.	+	+	+	+	+	+	+	+	+	+	+	+	+	-	+	+	+	+	+	+	+	+

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

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