



# MIC Test Strip Technical Sheet *Haemophilus influenzae*

## Specimen

Blood, CSF, sterile site (joint fluid, eye) and respiratory (sputum, tracheal aspirate, middle ear, nasopharynx)

<b>Medium</b>	Haemophilus Test Medium, Ref. 10080 or Mueller Hinton Fastidious Agar (Horse blood 5% + 20 mg/L $\beta$ -NAD), Ref. 10132
<b>Inoculum</b>	Suspension in broth to 0.5 McFarland (Ref. 80400) (if mucoid: 1 McFarland, Ref. 80401)
<b>Incubation</b>	35 $\pm$ 2 °C / 5% CO <sub>2</sub> / 20-24 hours
<b>Evaluating the results</b>	Bactericidal drugs: interpret the M.I.C. at complete inhibition of all growth including microcolonies, hazes and isolated colonies. Bacteriostatic drugs: interpret the M.I.C. at 80% inhibition when trailing is seen.

	Quality Control (M.I.C. $\mu$ g/mL)		CLSI INTERPRETATION M.I.C. Criteria ( $\mu$ g/mL)			EUCAST INTERPRETATION M.I.C. Criteria ( $\mu$ g/mL)		Example of ANTIBIOGRAM 140 mm petri dish
	<i>H. influenzae</i> ATCC® 49247	<i>H. influenzae</i> ATCC® 49766	S	I	R	S	R	
<b>AMP</b> AMPICILLIN	2-8	-	$\leq$ 1	2	$\geq$ 4	$\leq$ 1	>1	
<b>AUG</b> AMOXICILLIN/CLAVULANIC ACID 2/1 <sup>1</sup>	2-16	-	$\leq$ 4	-	$\geq$ 8			✓ or AMP
<b>AMC</b> AMOXICILLIN/CLAVULANIC ACID 2 $\mu$ g/mL <sup>1</sup>	-	-				$\leq$ 2	>2	
<b>AZM</b> AZITHROMYCIN (Ambient)	1-4	-	$\leq$ 4	-	-	$\leq$ 0.12	>4	
<b>AZM</b> AZITHROMYCIN (CO <sub>2</sub> ) <sup>2</sup>	4-16	-	$\leq$ 8	-	-			
<b>CTX</b> CEFOTAXIME	0.12-0.5	-	$\leq$ 2	-	-	$\leq$ 0.12	>0.12	✓ or CRO
<b>CRO</b> CEFTRIAXONE	0.06-0.25	-	$\leq$ 2	-	-	$\leq$ 0.12	>0.12	
<b>CXM</b> CEFUROXIME	-	0.25-1	$\leq$ 4	8	$\geq$ 16	$\leq$ 1	>2	
<b>CXM</b> CEFUROXIME oral	-	0.25-1				$\leq$ 0.12	>1	
<b>C</b> CHLORAMPHENICOL	0.25-1	-	$\leq$ 2	4	$\geq$ 8	$\leq$ 2	>2	✓ or CXM or LEV
<b>CLR</b> CLARITHROMYCIN (Ambient)	4-16	-	$\leq$ 8	16	$\geq$ 32	$\leq$ 1	>32	
<b>CLR</b> CLARITHROMYCIN (CO <sub>2</sub> ) <sup>2</sup>	8-32	-	$\leq$ 16	32	$\geq$ 64			
<b>LEV</b> LEVOFLOXACIN	0.008-0.03	-	$\leq$ 2	-	-	$\leq$ 1	>1	
<b>MRP</b> MEROPENEM	-	0.03-0.12	$\leq$ 0.5	-	-			✓
<b>MRP</b> MEROPENEM (infections other than meningitis)	-	-				$\leq$ 2	>2	
<b>MRP</b> MEROPENEM (meningitis)	-	-				$\leq$ 0.25	>1	
<b>TE</b> TETRACYCLINE	4-32	-	$\leq$ 2	4	$\geq$ 8	$\leq$ 1	>2	✓ or AZM or CLR
<b>SXT</b> TRIMETHOPRIM/SULFAMETHOXAZOLE 1/19 <sup>1</sup>	0.03-0.25	-	$\leq$ 0.5	1-2	$\geq$ 4	$\leq$ 0.5	>1	✓

## Notes

- Value on the M.I.C. scale refers to the first component of the combination.
- CLSI broth microdilution uses ambient incubation while agar based methods use CO<sub>2</sub> incubation that causes a pH decrease and may affect activity of macrolides. Quality control ranges and interpretive criteria for M.I.C. test are adjusted for CO<sub>2</sub> incubation.

MIC Test Strip, Patent No. 1395483

**LIOFILCHEM® s.r.l.**  
Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net liofilchem@liofilchem.net

## References

- CLSI M100-S23, 2013. Performance Standards for Antimicrobial Susceptibility Testing.
- CLSI M7-A9, 2012. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically.
- EUCAST. Breakpoint tables for interpretation of MICs and zone diameters Version 3.0, January 2013.



MTS01  
Rev.3 / 13.03.2013



# MIC Test Strip Technical Sheet **Enterococci**

## Specimen

Blood, wounds and sterile sites (tissues, peritoneal fluid), urines (urosepsis)

	<b>Enterococci (general)</b>	<b>HLAR (High level Aminoglycoside Resistance)</b>	<b>VRE (Vancomycin Resistant Enterococci)</b>
<b>Medium</b>	Mueller Hinton II Agar (Ref. 10031)	Mueller Hinton II Agar (Ref. 10031) <b>MIC Test Strip range:</b> Gentamicin and Streptomycin 0.064-1024 µg/ mL	Brain Heart Infusion Agar (Ref. 10060)
<b>Inoculum</b>	Suspension in physiological solution to 0.5 McFarland turbidity (Ref. 80400)	Suspension in physiological solution to 0.5-1 McFarland (heavier inoculum preferable) (1 McFarland Ref. 80401)	Suspension in broth to 2 McFarland. Dispense 0.1 mL per 90 mm agar plate and streak evenly.
<b>Incubation</b>	35 °C, ambient, 16-20 hours; interpret vancomycin at 24 hours.	35 °C, ambient, 24 hours (interpret streptomycin at 48 hours).	35 °C, ambient, 24-48 hours; interpret at 48 hours.
<b>Evaluating the results</b>	Bactericidal drugs: interpret at complete inhibition of all growth including microcolonies, hazes and isolated colonies. Bacteriostatic drugs: interpret at 80% inhibition when trailing is seen.	Read at complete inhibition including microcolonies, hazes and isolated colonies.	Look for hazes, microcolonies and isolated colonies. Use a magnifying glass, oblique light and/or tilt the plate to look for all growth.

	Quality Control M.I.C. µg/mL  <i>E. faecalis</i> ATCC® 29212	CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)		Examples of ANTIBIOGRAM	
		S	I	R	S	R	Detection of VRE (BHI) 90 mm petri dish	Other antibiotics and HLAR (Mueller Hinton) Observe: QDA for <i>E. faecium</i> only. 140 mm petri dish
<b>AMP</b> AMPICILLIN	0.5-2	≤8	-	≥16	<4	>8		✓
<b>LEV</b> LEVOFLOXACIN	0.25-2	≤2	4	≥8				
<b>LNZ</b> LINEZOLID	1-4	≤2	4	≥8	<4	>4		✓
<b>TEC</b> TEICOPLANIN	0.25-1	≤8	16	≥32	<2	>2		
<b>TE</b> TETRACYCLINE	8-32	≤4	8	≥16				✓
<b>VA</b> VANCOMYCIN	1-4	≤4	8-16	≥32	<4	>4		
<b>HLAR</b>								
<b>CN</b> GENTAMICIN	4-16							✓
<b>S</b> STREPTOMYCIN	64-256							✓
<b>VRE (BHI)</b>								
<b>TEC</b> TEICOPLANIN	0.25-1						✓	
<b>VA</b> VANCOMYCIN	2-6						✓	

**Phenotype Interpretation**

<b>HLAR</b>	<b>Negative</b>	<b>Positive</b>
GENTAMICIN	≤512	>512
STREPTOMYCIN	≤512	>1024

<b>VRE (BHI)</b>				
<b>Phenotype</b>	<b>Vancomycin (µg/mL)</b>		<b>Teicoplanin (µg/mL)</b>	<b>Species</b>
VanA	≥32 (R)	and	≥16 (I-R)	<i>E. faecalis</i> <i>E. faecium</i>
VanB	≥8-256 (I-R)	and	≤4 (S)	<i>E. faecalis</i> <i>E. faecium</i>
VanC1	4-16 (S-I)	and	≤4 (S)	<i>E. gallinarum</i>
VanC2	4-16 (S-I)	and	≤4 (S)	<i>E. casseliflavus</i> <i>E. flavescens</i>
VanD	64 (R)	and	≤4 (S)	<i>E. faecium</i>
VanE	16 (I)	and	≤4 (S)	<i>E. faecalis</i>

**References**

- CLSI M100-S23, 2013. Performance Standards for Antimicrobial Susceptibility Testing.
- CLSI M7-A9, 2012. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically.
- EUCAST. Breakpoint tables for interpretation of MICs and zone diameters Version 3.0, January 2013.

MIC Test Strip, Patent No. 1395483

**LIOFILCHEM® s.r.l.**

Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
 Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net liofilchem@liofilchem.net



MTS05  
 Rev.2 / 13.03.2013



# MIC Test Strip Technical Sheet **Meningococci**

## Specimen

Oro-nasopharynges, blood and CSF.

<b>Medium</b>	Mueller Hinton II Agar (Sheep blood 5%) Ref. (10131)
<b>Inoculum</b>	Suspension in broth to 0.5 McFarland (Ref. 80400) from chocolate agar (inoculum from SBA will contain 50% fewer CFU/mL).
<b>Incubation</b>	35 ± 2 °C/ 5% CO <sub>2</sub> / 20-24 hours
<b>Reading precaution</b>	Bactericidal drugs: interpret the M.I.C. at the point of complete inhibition of growth (microcolonies, hazes and isolated colonies). Bacteriostatic drugs: interpret the M.I.C. of hazy zone edges at 80% inhibition.

		Quality Control (M.I.C. µg/mL)		CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)		Example of ANTIBIOGRAM 140 mm petri dish
		<i>S. pneumoniae</i> ATCC® 49619	<i>E. coli</i> ATCC® 25922 <sup>1</sup>	S	I	R	S	R	
<b>P</b>	PENICILLIN G	0.25-1		≤0.06	0.12-0.25	≥0.5	≤0.06	>0.25	✓ 0.002-32
<b>CTX</b>	CEFOTAXIME	0.03-0.12		≤0.12			≤0.12	>0.12	
<b>CRO</b>	CEFTRIAZONE	0.03-0.12		≤0.12			≤0.12	>0.12	✓ or CTX
<b>CIP</b>	CIPROFLOXACIN		0.004-0.015	≤0.03	0.06	≥0.12	≤0.03	>0.06	✓ or MRP
<b>C</b>	CHLORAMPHENICOL	2-8		≤2	4	≥8	≤2	>4	
<b>MRP</b>	MEROPENEM	0.06-0.25		≤0.25			≤0.25	>0.25	
<b>RD</b>	RIFAMPICIN	0.015-0.06		≤0.5	1	≥2	≤0.25	>0.25	✓ 0.002-32
<b>SXT</b>	TRIMETHOPRIM/ SULFAMETHOXAZOLE 1/19 <sup>2</sup>	0.12-1		≤0.12	0.25	≥0.5	-	-	✓ or C

## Notes:

- CO<sub>2</sub> incubation.
  - Concentration refers to the first component of the combination.
- "<sup>1</sup>" indicates that susceptibility testing is not recommended as the species is a poor target for therapy with the drug.

## References


- CLSI M100-S23, 2013. Performance Standards for Antimicrobial Susceptibility Testing.
- CLSI M7-A9, 2012. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically.
- EUCAST. Breakpoint tables for interpretation of MICs and zone diameters Version 3.0 , January 2013.

MIC Test Strip, Patent No. 1395483



**LIOFILCHEM® s.r.l.**

Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net liofilchem@liofilchem.net



MTS08  
Rev.2 / 13.03.2013



# MIC Test Strip Technical Sheet Gram Negative Aerobes

Enterobacteriaceae, *Pseudomonas*, *Burkholderia*, *Acinetobacter* and *Stenotrophomonas* spp.

## Specimen

Blood, Cerebrospinal Fluid, sterile sites (joint, fluids, tissues), wounds, respiratory (sputum, transtracheal aspirate) and urines.

## Procedure

<b>Medium</b>	Mueller Hinton II Agar (ref. 10031).
<b>Inoculum</b>	Suspension in physiological solution to 0.5 McFarland (ref. 80400), mucoid strains: 1 McFarland (ref. 80401)
<b>Incubation</b>	35 ± 2 °C/ ambient / 16-20 hours non-fermenters: in case of low growth at 24 hours, confirm at 48 hours
<b>Evaluating the results</b>	Bactericidal drugs: interpret the M.I.C. at complete growth inhibition including microcolonies, hazes and isolated colonies. Bacteriostatic drugs: interpret the M.I.C. at 80% inhibition when trailing is seen.
<b>ESBL Extended Spectrum β-Lactamases</b>	Materials and procedure as above. Test Intensive Care Unit and critical isolates directly with MIC Test Strip CAZ/CAL and CTX/CTL strips. For other isolates, review aztreonam, cefotaxime, ceftazidime, ceftriaxone susceptibility results and use the interpretation criteria in the latest CLSI M100-S document. Confirm ESBL suspects using both MIC Test Strip CTX/CTL and CAZ/CAL. MIC Test Strip FEP/FEL may be used to test strains with non-determinable (ND) CTX/CTL and CAZ/CAL results.

		Quality Control (M.I.C. µg/mL)			CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)	
		<i>E. coli</i> ATCC® 25922	<i>P. aeruginosa</i> ATCC® 27853	<i>E. coli</i> ATCC® 35218	S	I	R	S	R
<b>AK</b>	AMIKACIN	0.5-4	1-4		≤16	32	≥64		
EUCAST	Enterobacteriaceae							≤8	>16
	<i>Pseudomonas</i> spp.							≤8	>16
	<i>Acinetobacter</i> spp.							≤8	>16
	Non-species related breakpoints							≤8	>16
<b>AMS</b>	AMPICILLIN/ SULBACTAM (2/1) <sup>1</sup>	2-8		8-32	≤8	16	≥32		
EUCAST	Enterobacteriaceae							≤8	>8
<b>ATM</b>	AZTREONAM	0.06-0.25	2-8						
CLSI	Enterobacteriaceae				≤4	8	≥16		
	<i>P. aeruginosa</i>				≤8	16	≥32		
	Non-Enterobacteriaceae				≤8	16	≥32		
EUCAST	Enterobacteriaceae							≤1	>4
	<i>Pseudomonas</i> spp.							≤1	>16
	Non-species related breakpoints							≤4	>8
<b>FEP</b>	CEFEPIME	0.015-0.12	0.5-4		≤8	16	≥32		
EUCAST	Enterobacteriaceae							≤1	>4
	<i>Pseudomonas</i> spp.							≤8	>8
	Non-species related breakpoints							≤4	>8

		Quality Control (M.I.C. µg/mL)			CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)									
		<i>E. coli</i> ATCC® 25922	<i>P. aeruginosa</i> ATCC® 27853	<i>E. coli</i> ATCC® 35218	S	I	R	S	R								
CTX	CEFOTAXIME	0.03-0.12	8-32														
	Enterobacteriaceae									≤1	2	≥4					
	CLSI <i>P. aeruginosa</i>									≤8	16-32	≥64					
	<i>Acinetobacter</i> spp.									≤8	16-32	≥64					
EUCAST	Non-Enterobacteriaceae				≤8	16-32	≥64										
	Enterobacteriaceae							≤1	>2								
	Non-species related breakpoints							≤1	>2								
CAZ	CEFTAZIDIME	0.06-0.5	1-4														
	Enterobacteriaceae									≤4	8	≥16					
	CLSI <i>P. aeruginosa</i>									≤8	16	≥32					
	<i>Acinetobacter</i> spp.									≤8	16	≥32					
	<i>B. cepacia</i>									≤8	16	≥32					
	<i>S. maltophilia</i>									≤8	16	≥32					
	Non-Enterobacteriaceae									≤8	16	≥32					
	Enterobacteriaceae															≤1	>4
	EUCAST <i>Pseudomonas</i> spp.															≤8	>8
	Non-species related breakpoints															≤4	>8
C	CHLORAMPHENICOL	2-8															
	EUCAST Enterobacteriaceae									≤8	16	≥32	≤8	>8			
CIP	CIPROFLOXACIN	0.004-0.015	0.25-1														
	Enterobacteriaceae (except <i>S. typhi</i> and extraintestinal <i>Salmonella</i> spp.)									≤1	2	≥4					
	CLSI <i>S. typhi</i> and extraintestinal <i>Salmonella</i> spp									≤0.06	0.12-0.5	1					
	<i>P. aeruginosa</i>									≤1	2	≥4					
	<i>Acinetobacter</i> spp.									≤1	2	≥4					
	Non-Enterobacteriaceae									≤1	2	≥4					
	Enterobacteriaceae														≤0.5	>1	
	EUCAST <i>Pseudomonas</i> spp.															≤0.5	>1
	<i>Acinetobacter</i> spp.															≤1	>1
	Non-species related breakpoints															≤0.5	>1
CS	COLISTIN	0.25-2	0.5-4														
	CLSI <i>P. aeruginosa</i>									≤2	4	≥8					
	<i>Acinetobacter</i> spp.									≤2	-	≥4					
	Non-Enterobacteriaceae									≤2	4	≥8					
	Enterobacteriaceae														≤2	>2	
	EUCAST <i>Pseudomonas</i> spp.															≤4	>4
<i>Acinetobacter</i> spp.							≤2	>2									

		Quality Control (M.I.C. µg/mL)			CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)	
		<i>E. coli</i> ATCC® 25922	<i>P. aeruginosa</i> ATCC® 27853	<i>E. coli</i> ATCC® 35218	S	I	R	S	R
<b>CN</b>	GENTAMICIN	0.25-1	0.5-2	≤4	8	≥16			
EUCAST	Enterobacteriaceae						≤2	>4	
	<i>Pseudomonas</i> spp.						≤4	>4	
	<i>Acinetobacter</i> spp.						≤4	>4	
	Non-species related breakpoints						≤2	>4	
<b>IMI</b>	IMIPENEM	0.06-0.25	1-4						
CLSI	Enterobacteriaceae			≤1	2	≥4			
	<i>P. aeruginosa</i>			≤2	4	≥8			
	<i>Acinetobacter</i> spp.			≤4	8	≥16			
	Non-Enterobacteriaceae			≤4	8	≥16			
EUCAST	Enterobacteriaceae						≤2	>8	
	<i>Pseudomonas</i> spp.						≤4	>8	
	<i>Acinetobacter</i> spp.						≤2	>8	
	Non-species related breakpoints						≤2	>8	
<b>LEV</b>	LEVOFLOXACIN	0.008-0.06	0.5-4						
CLSI	Enterobacteriaceae			≤0.12	0.25-1	≥2			
	<i>P. aeruginosa</i>			≤2	4	≥8			
	<i>Acinetobacter</i> spp.			≤2	4	≥8			
	<i>B. cepacia</i>			≤2	4	≥8			
	<i>S. maltophilia</i>			≤2	4	≥8			
	Non-Enterobacteriaceae			≤2	4	≥8			
EUCAST	Enterobacteriaceae						≤1	>2	
	<i>Pseudomonas</i> spp.						≤1	>2	
	<i>Acinetobacter</i> spp.						≤1	>2	
	Non-species related breakpoints						≤1	>2	
<b>MRP</b>	MEROPENEM	0.008-0.06	0.25-1						
CLSI	Enterobacteriaceae			≤1	2	≥4			
	<i>P. aeruginosa</i>			≤2	4	≥8			
	<i>Acinetobacter</i> spp.			≤4	8	≥16			
	Non-Enterobacteriaceae			≤4	8	≥16			
EUCAST	Enterobacteriaceae						≤2	>8	
	<i>Pseudomonas</i> spp.						≤2	>8	
	<i>Acinetobacter</i> spp.						≤2	>8	
	Non-species related breakpoints						≤2	>8	

		Quality Control (M.I.C. µg/mL)			CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)	
		<i>E. coli</i> ATCC® 25922	<i>P. aeruginosa</i> ATCC® 27853	<i>E. coli</i> ATCC® 35218	S	I	R	S	R
<b>TZP</b>	PIPERACILLIN / TAZOBACTAM (4 µg/mL) <sup>1</sup>	1-4	1-8	0.5-2					
CLSI	Enterobacteriaceae				≤16	32-64	≥128		
	<i>P. aeruginosa</i>				≤64	-	≥128		
	<i>Acinetobacter</i> spp.				≤16	32-64	≥128		
	Non-Enterobacteriaceae				≤16	32-64	≥128		
EUCAST	Enterobacteriaceae							≤8	>16
	<i>Pseudomonas</i> spp.							≤16	>16
	Non-species related breakpoints							≤4	>16
<b>PB</b>	POLYMYXIN B	0.25-2	1-4						
CLSI	<i>P. aeruginosa</i>				≤2	4	≥8		
	<i>Acinetobacter</i> spp.				≤2	-	≥4		
	Non-Enterobacteriaceae				≤2	4	≥8		
<b>TE</b>	TETRACYCLINE	0.5-2	8-32						
					≤4	8	≥16		
<b>TTC</b>	TICARCILLIN / CLAVULANIC ACID (2 µg/mL) <sup>1</sup>	4-16	8-32	8-32					
CLSI	Enterobacteriaceae				≤16	32-64	≥128		
	<i>P. aeruginosa</i>				≤64	-	≥128		
	<i>Acinetobacter</i> spp.				≤16	32-64	≥128		
	<i>B. cepacia</i>				≤16	32-64	≥128		
	<i>S. maltophilia</i>				≤16	32-64	≥128		
	Non-Enterobacteriaceae				≤16	32-64	≥128		
EUCAST	Enterobacteriaceae							≤8	>16
	<i>Pseudomonas</i> spp.							≤16	>16
	Non-species related breakpoints							≤8	>16
<b>SXT</b>	TRIMETHOPRIM / SULFAMETHOXAZOLE (1/19) <sup>1</sup>	≤0.5	8-32						
EUCAST	Enterobacteriaceae				≤2	-	≥4		>4
	<i>Pseudomonas</i> spp.							≤4	>4
	<i>Acinetobacter</i> spp.							≤2	>4

**Notes:**

1. Value on the M.I.C. scale refers to the first component of the combination.



		Examples of ANTIBIOGRAM					Muroid organisms e.g. <i>Klebsiella</i> spp., <i>Enterobacter</i> spp. and <i>P. aeruginosa</i>
		Enterobacteriaceae	<i>Acinetobacter</i> spp.	<i>Burkholderia</i> spp.	<i>Pseudomonas</i> spp.	<i>Stenotrophomonas</i> spp.	
		140 mm petri dish	140 mm petri dish	140 mm petri dish	140 mm petri dish	140 mm petri dish	140 mm petri dish
<b>AK</b>	AMIKACIN		✓ or CN		✓ or CN	✓	✓
<b>AMS</b>	AMPICILLIN/ SULBACTAM (2/1)		✓ or TTC			✓ or TTC	
<b>ATM</b>	AZTREONAM	✓			✓		
<b>C</b>	CHLORAMPHENICOL						
<b>CIP</b>	CIPROFLOXACIN	✓ or LEV	✓ or LEV		✓ or LEV	✓ or LEV	✓ or LEV
<b>CS</b>	COLISTIN						
<b>CN</b>	GENTAMICIN	✓ 0.016 - 256 or AK					
<b>IMI</b>	IMIPENEM	✓	✓ or MRP		✓ or MRP		✓
<b>LEV</b>	LEVOFLOXACIN			✓			
<b>MRP</b>	MEROPENEM			✓			
<b>TZP</b>	PIPERACILLIN / TAZOBACTAM (4 µg/mL)	✓			✓		✓
<b>TTC</b>	TICARCILLIN / CLAVULANIC ACID (2 µg/mL)						
<b>SXT</b>	TRIMETHOPRIM / SULFAMETHOXAZOLE (1/19)			✓		✓	
<b>CAZ</b>	CEFTAZIDIME			✓			✓
<b>CTX</b>	CEFOTAXIME						
<b>FEP</b>	CEFEPIME	✓ or CTX	✓ or CAZ		✓ or CAZ	✓ or CAZ	✓

		Quality Control (M.I.C. µg/mL)		Examples of ANTIBIOGRAM	
		<i>E. coli</i> ATCC® 35218 ESBL negative QC strain	<i>K. pneumoniae</i> ATCC® 700603 ESBL positive QC strain	Detection of ESBL 90 mm petri dish	Strains with ND results by CTX/CTL and CAZ/CAL 90 mm petri dish
<b>CAZ</b>	CEFTAZIDIME <sup>2</sup>	≤ 0.5	8-≥32		
<b>CAL</b>	CEFTAZIDIME + CLAV. ACID <sup>2</sup>	≤ 0.125	0.125-1		
<b>CTX</b>	CEFOTAXIME <sup>2,3</sup>	≤ 0.25	0.5-2		
<b>CTL</b>	CEFOTAXIME + CLAV. ACID <sup>2</sup>	0.016-0.064	0.125-1		
<b>FEP</b>	CEFEPIME <sup>2,3</sup>	≤ 0.25	0.5-2		
<b>FEL</b>	CEFEPIME + CLAV. ACID <sup>2</sup>	≤ 0.064	0.064-0.25		
<b>CTX/CTL</b>	CEFOTAXIME / CEFOTAXIME + CLAVULANIC ACID (4 µg/mL)			✓	
<b>CAZ/CAL</b>	CEFTAZIDIME / CEFTAZIDIME + CLAVULANIC ACID (4 µg/mL)			✓	
<b>FEP/FEL</b>	CEFEPIME / CEFEPIME + CLAVULANIC ACID (4 µg/mL)				✓

**Notes:**

- M.I.C. value below the strip range.
- Deformation of the ellipse is indicative of ESBL production even if the CTX/CTL or FEP/FEL ratio is <8.

**ESBL Phenotype Interpretation**

Negative	M.I.C. ratio of both CAZ/CAL and CTX/CTL	<8
Positive	M.I.C. for CTX $\geq 0.5$ and CTX/CTL ratio	$\geq 8$ OR
Positive	M.I.C. for CAZ $\geq 1$ and CAZ/CAL ratio	$\geq 8$ OR
Positive	M.I.C. for FEP/FEL	$\geq 8$
<b>IMPORTANT</b> Positive	"Phantom" zone or distortion of the CTX, CAZ or FEP inhibition ellipse confirms ESBL production, even if the CAZ/CAL, CTX/CTL or FEP/FEL ratio is < 8.	
Non-determinable (ND)	Off-scale results for both CTX/CTL and CAZ/CAL or one negative and the other off-scale. Strains with ND results for CTX/CTL and CAZ/CAL may be tested with FEP/FEL.	

**References**

- CLSI M100-S23, 2013. Performance Standards for Antimicrobial Susceptibility Testing.
- CLSI M7-A9, 2012. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically.
- EUCAST. Breakpoint tables for interpretation of MICs and zone diameters Version 3.0, January 2013.
- Rossolini, G.M. et al. (2011). Evaluation of a new gradient-diffusion system for MIC determination with Gram-negative pathogens. ECCMID, poster 572.

MIC Test Strip, Patent No. 1395483

**LIOFILCHEM® s.r.l.**

Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
 Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net liofilchem@liofilchem.net



MTS09  
 Rev.3 / 15.03.2013



# MIC Test Strip Technical Sheet **Gonococci**

## Specimen

Uro-genital tract, oropharynx, conjunctiva, blood, CSF and sterile sites and tissues.

<b>Medium</b>	GC agar base + defined supplements (CLSI) or Mueller Hinton Chocolate Agar (Ref. 10355)
<b>Inoculum</b>	Suspension in broth to 0.5 McFarland (Ref. 80400)
<b>Incubation</b>	36 ± 1 °C / 5% CO <sub>2</sub> / 20-24 hours
<b>Reading precautions</b>	Bactericidal drugs: interpret the M.I.C. at complete inhibition of all growth including microcolonies, hazes and isolated colonies. Bacteriostatic drugs: interpret the M.I.C. at 80% inhibition when trailing is seen.

		Quality Control (M.I.C. µg/mL) <i>N. gonorrhoeae</i> ATCC® 49226	CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)		Example of ANTIBIOGRAM 140 mm petri dish
			S	I	R	S	R	
<b>P</b>	PENICILLIN G <sup>1</sup>	0.25-1	≤0.06	0.12-1	≥2	≤0.06	>1	✓
<b>CRO</b>	CEFTRIAZONE	0.004-0.015	≤0.25			≤0.12	>0.12	✓
<b>CIP</b>	CIPROFLOXACIN	0.001-0.008	≤0.06	0.12-0.5	≥1	≤0.03	>0.06	✓
<b>TE</b>	TETRACYCLINE	0.25-1	≤0.25	0.5-1	≥2	≤0.5	>1	✓

## Notes:

1. If beta-lactamase positive, the strain should be reported as resistant to penicillins.

## References

- CLSI M100-S23, 2013. Performance Standards for Antimicrobial Susceptibility Testing.
- CLSI M7-A9, 2012. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically.
- EUCAST. Breakpoint tables for interpretation of MICs and zone diameters, Version 3.0, January 2013.

MIC Test Strip, Patent No. 1395483



**LIOFILCHEM® s.r.l.**

Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net liofilchem@liofilchem.net



MTS13  
Rev.2 / 12.02.2013



# MIC Test Strip Technical Sheet **Anaerobes**

## Specimen

Blood, wounds, respiratory (transtracheal aspirate) and sterile sites (CNS, tissues, joint fluids).

<b>Medium</b>	Brucella Blood Agar w Hemin and Vitamin K1 (ref. 10245)
<b>Inoculum</b>	Suspension in Brucella or "anaerobic" broth to 1 McFarland (Ref. 80401) turbidity. Make sure anaerobic conditions are kept at all times for obligate anaerobes.
<b>Incubation</b>	35 °C / anaerobic system / 24-72 hours. Obtain rapid anaerobiosis within 1-2 hours for Metronidazole.
<b>Evaluating the results</b>	Bactericidal drugs: interpret the M.I.C. at complete inhibition of growth including microcolonies, hazes and isolated colonies. Bacteriostatic drugs: interpret the M.I.C. at 80% inhibition when trailing is seen. If a "dip" is observed with MIC Test Strip Clindamycin, extrapolate the ellipse to the strip.

		Quality Control (M.I.C. µg/mL)			CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)		Example of ANTIBIOGRAM  140 mm petri dish
		<i>B. fragilis</i> ATCC® 25285	<i>B. thetaiotaomicron</i> ATCC® 29741	<i>E. lentum</i> ATCC® 43055	S	I	R	S	R	
<b>AUG</b>	AMOXICILLIN/CLAVULANIC ACID 2/1 <sup>1)</sup>	0.25-1	0.5-2	-	≤ 4	8	≥ 16	≤ 4	≥ 8	
<b>AMS</b>	AMPICILLIN/ SULBACTAM <sup>1)</sup>	0.5-2	0.5-2	0.5-2	≤ 8	16	≥ 32	≤ 4	≥ 8	
<b>P</b>	PENICILLIN G	8-32	8-32	-	≤ 0.5	1	≥ 2	≤ 0.25	> 0.5	✓
<b>FOX</b>	CEFOXITIN	2-8	8-64	2-16	≤ 16	32	≥ 64			✓
<b>CD</b>	CLINDAMYCIN	0.5-2	2-8	0.06-0.25	≤ 2	4	≥ 8	≤ 4	> 4	✓
<b>ETP</b>	ERTAPENEM	0.06-0.25	0.5-2	0.5-4	≤ 4	8	≥ 16	≤ 1	> 1	✓ or MRP or IMI
<b>IMI</b>	IMIPENEM	0.03-0.12	0.25-1	0.25-2	≤ 4	8	≥ 16	≤ 2	> 8	
<b>MRP</b>	MEROPENEM	0.03-0.25	0.06-0.5	0.125-1	≤ 4	8	≥ 16	≤ 2	> 8	
<b>LZ</b>	METRONIDAZOLE <i>C. difficile</i>	0.25-1	0.5-2	-	≤ 8	16	≥ 32	≤ 4 ≤ 2	> 4 > 2	✓
<b>TZP</b>	PIPERACILLIN/TAZOBACTAM 4 <sup>1)</sup>	0.03-0.25	2-16	8-32	≤ 32	64	≥ 128	≤ 8	> 16	✓ or AUG or AMS or TTC
<b>TTC</b>	TICARCILLIN/CLAVULANIC ACID 2 <sup>1)</sup>	0.06-0.5	0.5-2	8-32	≤ 32	64	≥ 128	≤ 8	> 16	

## Notes:

1. Value on MIC scale refers to the first component of the combination.

## References

- CLSI M11-A7. Methods for Dilution Antimicrobial Susceptibility Testing of Anaerobic Bacteria.
- CLSI M100-S23. Performance Standards for Antimicrobial Susceptibility Testing, January 2013.
- EUCAST. Breakpoint tables for interpretation of MICs and zone diameters, Version 3.1, February 2013.

MIC Test Strip, Patent No. 1395483

**LIOFILCHEM® s.r.l.**  
Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net liofilchem@liofilchem.net



MTS14  
Rev.4 / 29.05.2013



# MIC Test Strip Technical Sheet *Campylobacter* spp.

## Specimen

Stools, blood, tissues.

<b>Medium</b>	Mueller Hinton II Agar (Sheep blood 5%) (Ref. 10131)
<b>Inoculum</b>	Suspension of 24-48 hours colonies from blood agar in broth to 1 McFarland (Ref. 80401)
<b>Incubation</b>	Microaerophilic incubation, 10% CO <sub>2</sub> / 5% O <sub>2</sub> / 85% N <sub>2</sub> Do not invert plates due to excessive capsular material that may be produced. <i>C. jejuni</i> and <i>C. coli</i> : 42 °C for 24 hours; 35 °C for 48-72 hours. <i>C. fetus</i> 35 °C for 48 hours.
<b>Evaluating the results</b>	No CLSI criteria available. Use publications for suggestions of possible interpretive criteria.
<b>Reading Precaution</b>	<i>Campylobacter</i> colonies may be translucent and difficult to interpret. Tilt the plate and/or use oblique light or a magnifying glass when reading the M.I.C. end point. Capsular material from highly mucoid strains may deform the inhibition ellipse. Repeat the test if necessary and do not incubate plates upside down. Bactericidal drugs: interpret at complete inhibition of all growth, including microcolonies, hazes and isolated colonies. Bacteriostatic drugs: interpret at 80% inhibition when trailing is seen.

	Quality Control M.I.C. µg/mL		EUCAST INTERPRETATION M.I.C. µg/mL		Example of ANTIBIOGRAM  140 mm petri dish
	<i>C. jejuni</i> ATCC® 33560 35 °C for 48 hours	<i>C. jejuni</i> ATCC® 33560 42 °C for 24 hours <sup>1</sup>	S	R	
<b>AZM</b> AZITHROMYCIN	0.03-0.25	0.03-0.12			
<b>CIP</b> CIPROFLOXACIN	0.06-0.25	0.03-0.12	≤0.5	>0.5	✓
<b>DX</b> DOXYCYCLINE	0.12-0.5	0.12-0.5			✓ or MRP
<b>E</b> ERYTHROMYCIN	0.5-2	0.25-2	≤4	>4	✓
<b>CN</b> GENTAMICIN	0.5-2	0.25-2			✓ 0.016 - 256
<b>LEV</b> LEVOFLOXACIN	0.06-0.25	0.03-0.25			
<b>MRP</b> MEROPENEM	0.008-0.03	0.008-0.03			
<b>TE</b> TETRACYCLINE	0.25-2	0.25-2	≤2	>2	

## Notes


1. Some clinical isolates of *C. jejuni* spp. *doyley*, *C. fetus* and *C. lari* may not grow at 42 °C. Test these isolates at 35 °C.

## References

- CLSI M45-A, 2006. Methods for Antimicrobial Dilutions and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline.
- EUCAST. Breakpoint tables for interpretation of MICs and zone diameters, Version 3.0, January 2013.

MIC Test Strip, Patent No. 1395483

**LIOFILCHEM® s.r.l.**  
Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net liofilchem@liofilchem.net



MTS16  
Rev.2 / 15.03.2013



# MIC Test Strip Technical Sheet **Streptococci**

## Specimen

Blood, CSF, wounds, sterile sites (joint fluid, eye, tissues) and respiratory samples (sputum, tracheal aspirate, middle ear fluid, nasopharynx).

	<b>Streptococci</b>	<b>Abiotrophia and Granulicatella spp.</b>
<b>Medium</b>	Mueller Hinton II Agar (Sheep blood 5%), Ref. 10131 or Mueller Hinton Fastidious Agar (Horse blood 5% + 20 mg/L β-NAD), Ref. 10132	Mueller Hinton Chocolate + 0.001% pyridoxal HCl + 0.01% cysteine IsoSensitest + 5% human blood + 0.001% pyridoxal HCl + 0.01% cysteine
<b>Inoculum</b>	Suspension in broth to 0.5 McFarland (Ref. 80400)	Suspension in broth to 1 McFarland (Ref. 80401)
<b>Incubation</b>	35 ± 2 °C / 5% CO <sub>2</sub> / 20- 24 hours	35 ± 2 °C / 5% CO <sub>2</sub> / 20- 24 hours
<b>Interpretation</b>	Bactericidal drugs: interpret the M.I.C. at complete inhibition of growth including microcolonies, hazes and isolated colonies. Bacteriostatic drugs: interpret the M.I.C. at 80% inhibition when trailing is seen.	

		<b>Quality Control</b> (M.I.C. µg/mL) <i>S. pneumoniae</i> ATCC® 49619	<b>CLSI INTERPRETATION</b> M.I.C. Criteria (µg/mL)			<b>EUCAST INTERPRETATION</b> <sup>1</sup> M.I.C. Criteria (µg/mL)		<b>Example of ANTIBIOGRAM</b> 140 mm petri dish
			<b>S</b>	<b>I</b>	<b>R</b>	<b>S</b>	<b>R</b>	
<b>P</b>	PENICILLIN G	0.25-1						✓ 0.002 - 32
CLSI	<i>Streptococcus</i> spp. β-Hemolytic Group		≤0.12	-	-			
	<i>Streptococcus</i> spp. Viridans		≤0.12	0.25-2	≥4			
EUCAST	Streptococcus groups A,B,C and G					≤0.25	>0.25	
	Viridans group Streptococci					≤0.25	>2	
<b>CTX</b>	CEFOTAXIME	0.03-0.12						✓ (or CRO)
CLSI	<i>Streptococcus</i> spp. β-Hemolytic Group		≤0.5	-	-			
	<i>Streptococcus</i> spp. Viridans		≤1	2	≥4			
EUCAST	Viridans group Streptococci					≤0.5	>0.5	
<b>CRO</b>	CEFTRIAXONE	0.03-0.12						
CLSI	<i>Streptococcus</i> spp. β-Hemolytic Group		≤0.5	-	-			
	<i>Streptococcus</i> spp. Viridans		≤1	2	≥4			
EUCAST	Viridans group Streptococci					≤0.5	>0.5	
<b>C</b>	CHLORAMPHENICOL	2-8	≤4	8	≥16	≤8	>8	✓ (or E)
<b>CD</b>	CLINDAMYCIN (ambient)	0.03-0.12	≤0.25	0.5	≥1	≤0.5	>0.5	✓ (or TE)
<b>CD</b>	CLINDAMYCIN (CO <sub>2</sub> )	0.064-0.25	≤0.5	1	≥2			
<b>DAP</b>	DAPTOMYCIN	0.064-0.5	≤1	-	-	≤1	>1	
<b>E</b>	ERYTHROMYCIN (ambient)	0.032-0.125	≤0.25	0.5	≥1	≤0.25	>0.5	
<b>E</b>	ERYTHROMYCIN (CO <sub>2</sub> )	0.064-0.25	≤1	2	≥4			
<b>OFX</b>	OFLOXACIN	1-4	≤2	4	≥8			✓
<b>TE</b>	TETRACYCLINE	0.06-0.5	≤2	4	≥8	≤1	>2	
<b>VA</b>	VANCOMICIN	0.12-0.5	≤1	-	-	≤2	>2	✓

**Note**


1. If not specified EUCAST interpretation refers to Streptococcus groups A,B,C and G.

**References**

- CLSI M100-S23, 2013. Performance Standards for Antimicrobial Susceptibility Testing.
- CLSI M7-A9, 2012. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically.
- EUCAST. Breakpoint tables for interpretation of MICs and zone diameters Version 3.0 , January 2013.

MIC Test Strip, Patent No. 1395483

 **LIOFILCHEM® s.r.l.**  
Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
Tel. +39 0858930745 Fax +39 0858930330 [www.liofilchem.net](http://www.liofilchem.net) [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)



MTS19  
Rev.3 / 15.03.2013



# MIC Test Strip Technical Sheet **Staphylococci**

## Specimen

Blood, wounds, sterile sites (tissues, bone, joints, fluids, CNS) and indwelling devices.

## Procedure

**Medium:** Mueller Hinton II Agar (ref. 10031).

**Inoculum:** Suspension in physiological solution to 0.5 McFarland (Ref. 80400).

**Incubation:** 35 ± 2 °C / ambient / 16-20 hours. Interpret vancomycin and oxacillin after 24 hours.

**Interpretation of results:** Bactericidal drugs: interpret the M.I.C. at complete inhibition of growth including microcolonies, hazes and isolated colonies. Bacteriostatic drugs: interpret the M.I.C. at 80% inhibition when trailing is seen.

**ORSA/ OR-CNS** (Oxacillin resistant *S. aureus* and Coagulase Negative Staphylococci).

**BORSA** (Bordeline Oxacillin Resistant *S. aureus*).

**Medium:** Mueller Hinton Agar + 2% NaCl (Ref. 11206).

**Inoculum:** Suspension in physiological solution to 0.5-1 McF (heavier inoculum improves detection of low level R).

**Incubation:** 35 ± 2 °C / ambient / 24 hours for ORSA/ BORSA, 48 hours for ORCNS.

**Interpretation of results:** Interpret at complete inhibition of all growth; watch for microcolonies, hazes and isolated colonies.

**GISA/hGISA** (Glycopeptide Intermediate / Resistant and Hetero-Intermediate / Resistant *S. aureus*).

**Medium:** Brain Heart Infusion Agar (Ref. 10160).

**Inoculum:** Suspension in broth to 2 McF (heavier inoculum improves detection of hetero-resistance).

**Incubation:** 35 ± 2 °C / ambient / interpret at 24 hours and confirm at 48 hours.

**Interpretation of results:** Interpret at complete inhibition; watch for hazes, microcolonies and isolated colonies. Use a magnifying glass, oblique light and tilt the plate.

		Quality Control (M.I.C. µg/mL)		CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)		Examples of ANTIBIOGRAM			
		<i>S. aureus</i> ATCC® 29213	<i>S. aureus</i> ATCC® 43300	S	I	R	S	R	ORSA/ORCNS/ BORSA (Mueller Hinton + 2% NaCl) 90 mm petri dish	Detection of glycopeptide resistance 90 mm petri dish	For confirmed ORSA/ORCNS (Muller Hinton Agar) 140 mm petri dish	For Non-ORSA/ ORCNS (Mueller Hinton Agar) 140 mm petri dish
<b>AUG</b>	AMOXICILLIN/CLAVULANIC ACID 2/1 <sup>1</sup>	0.12-0.5		≤4	-	≥8			✓			
<b>P</b>	PENICILLIN G	0.25-2		≤0.12	-	≥0.25	≤0.12	>0.12				✓
<b>C</b>	CHLORAMPHENICOL	2-16		≤8	16	≥32	≤8	>8				
<b>CIP</b>	CIPROFLOXACIN	0.12-0.5		≤1	2	≥4	≤1	>1				✓
<b>CD</b>	CLINDAMYCIN	0.06-0.25		≤0.5	1-2	≥4	≤0.25	>0.5			✓	✓
<b>DAP</b>	DAPTOMYCIN	0.12-1		≤1	-	-	≤1	>1			✓	
<b>E</b>	ERYTHROMYCIN	0.25-1		≤0.5	1-4	≥8	≤1	>2				✓
<b>CN</b>	GENTAMICIN	0.12-1		≤4	8	≥16	≤1	>1				
<b>LNZ</b>	LINEZOLID	1-4		≤4	-	≥8	≤4	>4			✓	✓
<b>RD</b>	RIFAMPICIN	0.004-0.015		≤1	2	≥4	≤0.06	>0.5				
<b>TEC</b>	TEICOPLANIN <i>S. aureus</i> Coagulase Negative Staphylococci	0.25-1		≤8	16	≥32	≤2 ≤4	>2 >4				
<b>TE</b>	TETRACYCLINE	0.12-1		≤4	8	≥16	≤1	>2				
<b>TGC</b>	TIGECYCLINE	0.03-0.25					≤0.5	>0.5				



	Quality Control (M.I.C. µg/mL)	CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)		Examples of ANTIBIOGRAM					
		<i>S. aureus</i> ATCC® 29213	<i>S. aureus</i> ATCC® 43300	S	I	R	S	R	ORSA/ORCNS/ BORSA (Mueller Hinton + 2% NaCl) 90 mm petri dish	Detection of glycopeptide resistance 90 mm petri dish	For confirmed ORSA/ORCNS (Muller Hinton Agar) 140 mm petri dish	For Non-ORSA/ ORCNS (Mueller Hinton Agar) 140 mm petri dish
<b>SXT</b> TRIMETHOPRIM / SULFAMETHOXAZOLE 1/19 <sup>1</sup>	≤0.5			≤2	-	≥4	≤2	>4				✓
<b>VA</b> VANCOMICIN <i>S. aureus</i> <i>Staphylococcus</i> spp. Coagulase Negative Staphylococci	0.5-2			≤2 ≤4	4-8 8-16	≥16 ≥32	≤2 ≤4	>2 >4			✓	
<b>ORSA</b>												
<b>OX</b> OXACILLIN <i>S. aureus</i> and <i>S. lugdunensis</i> Coagulase Negative Staphylococci except <i>S. lugdunensis</i>	0.12-0.5	16 - 24		≤2 ≤0.25	- -	≥4 ≥0.5			✓			
<b>GISA/hGISA</b>												
<b>TEC</b> TEICOPLANIN	0.5-2									✓		
<b>VA</b> VANCOMICIN	1-4									✓		

#### Phenotype Interpretation

	OXACILLIN	AMOXICILLIN/CLAVULANIC ACID
<b>ORSA (mecA+)</b>	≥4	≥8 <sup>1</sup>
<b>ORCNS (mecA+)</b>	≥0.5	≥8 <sup>1</sup>
<b>BORSA (mecA-)</b>	≥4	≤4 <sup>1</sup>
	VANCOMICIN	TEICOPLANIN
<b>GISA/hGISA</b>	≥8	AND ≥8 OR ≥12

#### Notes

1. Value on the M.I.C. scale refers to the first component of the combination.

#### References

- CLSI M100-S23, 2013. Performance Standards for Antimicrobial Susceptibility Testing.
- CLSI M7-A9, 2012. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically.
- EUCAST. Breakpoint tables for interpretation of MICs and zone diameters Version 3.0, January 2013.
- Stefani, S. et al (2011). A new reliable screening method for the evaluation of VISA and hVISA strains by "Vancomycin-Teicoplanin MIC Test Strip" (VTMTS). ECCMID poster 776.

MIC Test Strip, Patent No. 1395483

**LIOFILCHEM® s.r.l.**  
Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net liofilchem@liofilchem.net



MTS20  
Rev.2 / 15.03.2013



# MIC Test Strip Technical Sheet *Helicobacter pylori*

## Specimen

Stomach biopsy, composed of samples from multiple sites.

<b>Medium</b>	Mueller Hinton II Agar (Sheep blood 5%) (Ref. 10131)
<b>Inoculum</b>	72 h (or older) viable colonies are suspended in broth (Mueller Hinton or other) supplemented with 5% serum; adjust turbidity to 3 McFarland. Use 1 MIC Test Strip per 90 mm plate; position the handle of the strip against the edge of the plate.
<b>Incubation</b>	35 ± 2 °C/ microaerophilic (atmosphere produced by a gas-generating system suitable for Campylobacter) 72 hours (or longer i.e. until a visible inhibition ellipse is seen). For metronidazole, a 24 hours anaerobic pre-incubation followed by 48 hours or longer microaerophilic incubation has been recommended by some investigators as a better option.
<b>Reading precautions</b>	<i>H. pylori</i> colonies are pin-point, translucent and difficult to see. Tilt the plate and/or use oblique light or a magnifying glass when reading the M.I.C. endpoint. Bactericidal drugs: interpret the M.I.C. at the point of complete inhibition of all growth, including microcolonies, hazes and isolated colonies. Bacteriostatic drugs: interpret the M.I.C. of hazy zone edges at 80% inhibition.

	Quality Control (M.I.C. µg/mL) <i>H. pylori</i> ATCC® 43504	CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)		Example of ANTIBIOGRAM 140 mm petri dish
		S	I	R	S	R	
<b>AML</b> AMOXICILLIN	0.015-0.12	≤1	-	≥1	≤0.12	>0.12	✓
<b>CLR</b> CLARITHROMYCIN	0.015-0.12	≤0.25	5	≥1	≤0.25	>0.5	✓
<b>LZ</b> METRONIDAZOLE	64-256	≤4	-	≥4	≤8	>8	✓
<b>TE</b> TETRACYCLINE	0.12-1	≤2	-	≥2	≤1	>1	✓

## References

- CLSI M100-S20, 2010. Performance Standards for Antimicrobial Susceptibility Testing.
- EUCAST. Breakpoint tables for interpretation of MICs and zone diameters Version 3.0, January 2013.

MIC Test Strip, Patent No. 1395483

**LIOFILCHEM® s.r.l.**  
Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net liofilchem@liofilchem.net



MTS21  
Rev.2 / 15.03.2013



## MIC Test Strip Technical Sheet *Streptococcus pneumoniae*

### Specimen

Blood, CSF and respiratory sites (sputum, tracheal aspirate, middle ear fluid, nasopharynx).

<b>Medium</b>	Mueller Hinton II Agar (Sheep blood 5%), Ref. 10131 or Mueller Hinton Fastidious Agar (Horse blood 5% + 20 mg/L β-NAD), Ref. 10132.
<b>Inoculum</b>	Suspension in Mueller Hinton Broth (Ref. 24107) to 0.5 McFarland (Ref. 80400). 1 McFarland (Ref. 80401) for mucoid.
<b>Incubation</b>	35 ± 2 °C/ 5% CO <sub>2</sub> / 20-24 hours.
<b>Interpretation of results</b>	Bactericidal drugs: interpret the M.I.C. at complete inhibition of growth including microcolonies, hazes and isolated colonies. Bacteriostatic drugs: interpret the M.I.C. at 80% inhibition when trailing is seen.

		Quality Control (M.I.C. µg/mL) <i>S. pneumoniae</i> ATCC® 49619	CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)		Examples of ANTIBIOGRAM	
			S	I	R	S	R	140 mm petri dish	90 mm petri dish
<b>AZM</b>	AZITHROMYCIN (Ambient)	0.06-0.25	≤0.5	1	≥2	≤0.25	>0.5		
<b>AZM</b>	AZITHROMYCIN (CO <sub>2</sub> ) <sup>1</sup>	0.5-2	≤4	8	≥16				
<b>P</b>	PENICILLIN G	0.25-1						✓ 0.002 - 32	✓ 0.002 - 32
	Parenteral (non-meningitis)		≤2	4	≥8	≤0.06	>2		
	Parenteral (meningitis)		≤0.06	-	≥0.12	≤0.06	>0.06		
	Oral		≤0.06	0.12-1	≥2				
<b>CTX</b>	CEFOTAXIME	0.03-0.12							
	meningitis		≤0.5	1	≥2	≤0.5	>2	✓ (or CRO)	✓ (or CRO)
	non-meningitis		≤1	2	≥4				
<b>CRO</b>	CEFTRIAXONE	0.03-0.12							
	meningitis		≤0.5	1	≥2	≤0.5	>2	✓ (or CTX)	✓ (or CTX)
	non-meningitis		≤1	2	≥4				
<b>C</b>	CHLORAMPHENICOL	2-8	≤4	-	≥8	≤8	>8		
<b>CLR</b>	CLARITHROMYCIN (Ambient)	0.03-0.12	≤0.25	0.5	≥1	≤0.25	>0.5		
<b>CLR</b>	CLARITHROMYCIN (CO <sub>2</sub> ) <sup>1</sup>	0.064-0.25	≤0.5	1	≥2				
<b>CD</b>	CLINDAMYCIN (Ambient)	0.03-0.12	≤0.25	0.5	≥1	≤0.5	>0.5		
<b>CD</b>	CLINDAMYCIN (CO <sub>2</sub> ) <sup>1</sup>	0.064-0.25	≤0.5	1	≥2				
<b>ETP</b>	ERTAPENEM	0.03-0.25	≤1	2	≥4	≤0.5	>0.5		
<b>E</b>	ERYTHROMYCIN (Ambient)	0.03-0.12	≤0.25	0.5	≥1	≤0.25	>0.5	✓	
<b>E</b>	ERYTHROMYCIN (CO <sub>2</sub> ) <sup>1</sup>	0.064-0.25	≤1	2	≥4				
<b>IMI</b>	IMIPENEM	0.03-0.12	≤0.12	0.25-0.5	≥1	≤2	>2		
<b>LEV</b>	LEVOFLOXACIN	0.5-2	≤2	4	≥8	≤2	>2		
<b>LNZ</b>	LINEZOLID	0.25-2	≤2	-	-	≤2	>4		

	Quality Control (M.I.C. µg/mL) <i>S. pneumoniae</i> ATCC® 49619	CLSI INTERPRETATION M.I.C. Criteria (µg/mL)			EUCAST INTERPRETATION M.I.C. Criteria (µg/mL)		Examples of ANTIBIOGRAM	
		S	I	R	S	R	140 mm petri dish	90 mm petri dish
<b>MRP</b> MEROPENEM meningitis non-meningitis	0.06-0.25	≤0.25	0.5	≥1	≤0.25 ≤2	>1 >2	✓	
<b>OFX</b> OFLOXACIN	1-4	≤2	4	≥8	≤0.12	>4		
<b>TE</b> TETRACYCLINE	0.06-0.5	≤1	2	≥4	≤1	>2		
<b>SXT</b> TRIMETHOPRIM/ SULFAMETHOXAZOLE 1/19 <sup>2</sup>	0.12-1	≤0.5	1-2	≥4	≤1	>2	✓	
<b>VA</b> VANCOMYCIN	0.12-0.5	≤1	-	-	≤2	>2	✓	

**Notes**


1. CLSI broth microdilution uses ambient incubation while agar based methods use CO<sub>2</sub> incubation that causes a pH decrease and may affect activity of macrolides and lincosamides. Quality control ranges and interpretive criteria for MIC Test Strip are adjusted for CO<sub>2</sub> incubation.
2. Values on the MIC scale refer to the first component of the combination.

**References**

- CLSI M100-S23, 2013. Performance Standards for Antimicrobial Susceptibility Testing.
- CLSI M7-A9, 2012. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically.
- EUCAST. Breakpoint tables for interpretation of MICs and zone diameters Version 3.0, January 2013.

MIC Test Strip, Patent No. 1395483

 **LIOFILCHEM® s.r.l.**  
Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net liofilchem@liofilchem.net



MTS23  
Rev.4 / 15.03.2013



# MIC Test Strip Technical Sheet **Direct MIC**

## Direct MIC Testing of Critical Specimens

### INTENDED USE

Direct specimen testing may provide faster information for therapy guidance and/or correction of empiric therapy in urgent clinical situations. However, results with direct specimen testing must always be considered preliminary and the clinician cautioned until confirmed by standardised pure isolate testing.

### DIRECT SPECIMEN TESTING WITH MIC TEST STRIP

- Positive blood culture specimen testing with Gram positive and negative aerobes, anaerobes and yeast, have been investigated and published.
- Sputum testing for patients with cystic fibrosis and lower respiratory tract secretion samples from ventilator associated pneumonia (VAP) have also been studied.
- Variations in inoculum, type of organisms and potential contamination will not affect results significantly, as these phenomena can be visually inspected on the agar plate.

### TEST PROCEDURE

1. Specimens: CSF, urines, sputum, respiratory tract samples and positive blood cultures from critical infections and critical patients.
2. Perform Gram stain or India ink / Lacto-phenol cotton blue stain (yeast) and examine microscopically.
3. Use a rich media and different incubation conditions to cover different suspected organisms:

Gram positive aerobes:	Mueller Hinton Agar + 5% blood/ 35°C/ ambient and 5% CO <sub>2</sub> / 20-48h
Gram negative aerobes:	Mueller Hinton Agar + 5% blood/ 35°C/ ambient/ 20-48h
Anaerobes:	Plate 1 – Brucella agar + 5% blood + vitamin K (1 µg/mL) + hemin (5 µg/mL) (BBA) 35°C/ anaerobic system/ 24-48h  Plate 2 – Mueller Hinton agar +5% blood/ 35°C/ ambient/ 20-48h
Yeast:	RPMI / 35°C/ moist in plastic bag/ 24-48h

4. When testing sputum, sputulise the sample.
5. Pipette 0.3 mL of the undiluted positive blood culture, CSF or urine onto the agar plate and streak out evenly. If cells are sparse (microscopy), centrifuge to concentrate, re-suspend and streak. For sputum, moisten the swab with the sputolysed specimen and streak evenly onto the agar plate.
6. Test key Gram positive or Gram negative drugs as guided by Gram stain, microscopy and suspected/ expected pathogens. For yeast, test fluconazole, itraconazole, voriconazole and amphotericin.
7. Start examining agar plates for preliminary results if growth is clearly visible for rapidly growing aerobes after 6-8h, 12-16h, then confirm again at 24h or longer. For sputum, incubate 48-72h.
  - Resistant results are considered more useful.
  - Susceptible results should be treated with caution until confirmed by standardised testing.
8. For sputum, identify the different colony morphotypes, species, growth patterns and respective MTS endpoints. Document the interaction between pathogens and normal flora, take a photograph of the agar plate to document significant findings and discuss with the treating clinician.
9. Always inform and caution the clinician that:
  - A modified procedure was used to generate preliminary results only.
  - Final results from standardised testing are pending.
10. ALWAYS CONFIRM DIRECT TESTING RESULTS WITH THE STANDARDISED OVERNIGHT MTS PROCEDURE.
11. Initiate this procedure simultaneously with direct testing and report results as soon as available.
12. Collect data to validate rapid specimen versus pure isolate testing.
13. Run QC organisms on all test occasions.

<b>ANTIBIOGRAM – EXAMPLES ONLY (PLEASE USE YOUR OWN FORMULARY)</b>			
Gram positive aerobic diplococci	Penicillin G (P) 0.002 - 32 µg/mL Cefotaxime (CTX) 0.002 - 32 µg/mL	Meropenem (MRP) Vancomycin (VA)	Trimethoprim / Sulfamethoxazole (SXT) Erythromycin (E)
Gram positive aerobic cocci	Gatifloxacin (GAT) Cefoxitin (FOX)	Gentamicin (CN) Linezolid (LNZ)	Penicillin G (P) 0.016 - 256 µg/mL Vancomycin (VA)
Gram negative aerobic bacilli	Amikacin (AK) Ciprofloxacin (CIP)	Aztreonam (ATM) Imipenem (IMI)	Cefepime (FEP) Piperacillin / Tazobactam (TZP)
Anaerobes cocci/ bacilli	Penicillin G (P) 0.016 - 256 µg/mL Imipenem (IMI)	Cefoxitin (FOX) Metronidazole (LZ)	Clindamycin (CD) Piperacillin / Tazobactam (TZP)
Yeast	Fluconazole (FLU) Amphotericin B (AMB)	Itraconazole (ITC)	Voriconazole (VO)

#### REFERENCES

- CLSI M100-S22, 2012. Performance Standards for Antimicrobial Susceptibility Testing.
- CLSI M7-A9, 2012. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically.
- CLSI M11-A7, 2007. Methods for Dilution Antimicrobial Susceptibility Testing of Anaerobic Bacteria.
- CLSI M11-S1 Performance Standards for Antimicrobial Susceptibility Testing of Anaerobic Bacteria.
- CLSI M27-A3. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard - Third Edition.
- CLSI M27-S3. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Third Informational supplement.
- EUCAST. Breakpoint tables for interpretation of MICs and zone diameters Version 2.0, January 2012.
- Rossolini, G.M. et al. (2011). Evaluation of a new gradient-diffusion system for MIC determination with Gram-negative pathogens. ECCMID, poster 572.
- Stefani, S. et al (2011). A new reliable screening method for the evaluation of VISA and hVISA strains by "Vancomycin-Teicoplanin MIC Test Strip" (VTMTS). ECCMID poster 776.



**LIOFILCHEM® s.r.l.**

Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy  
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net liofilchem@liofilchem.net



MTS33  
Rev.0 / 29.03.2012