SAFE® READ Single Tube Mobile

Mobile 2D Barcode Single Tube Reader

- Intelligent and fully portable 2D barcode reader for single tubes
- Wired or wireless communications
- Up- or down-load data from host computer
- Operating modes selected from main screen icons
- Full day operation without the need to re-charge the battery
- Cryoprotection provided as standard
- Can be customised for applications in your laboratory
- Client software to provide wired interface between computer and instrument
- With Battery charger and USB lead to connect instrument to host computer
- Dimensions: 155mm x 65mm x 38mm
- Battery life: Typically, 8 hours between charges
- Data storage on instrument: Typically, up to 60.000 samples
- Small 1D linear code scanner linked to the Handheld by Bluetooth Protective Case
- Protective cover
- Optional for Wireless Applications: Handheld Server, Keyboard Wedge, Link to Excel, Link to Database

Item Number

DMST-H

Order Information



