

Technical Datasheet

Date: 02.01.25

Page 1 of 2

EN ISO 13485 certified		
Article - No:	06-365-2018	Revision: 11
Remarks		single use only

Specification

Product	Pipette tips PP, premium surface, 0-200µl, 96 pcs/rack, transparent,, moulded rings, np premium bio, sterile R, autocl. up to 121°C
Dimensions	length: 53,40 mm (+/- 0,2 mm)
	inner diameter: 0,47 mm (+/- 0,2 mm)
Colour	transparent,
Graduation	moulded rings
Nominal Capacity	0-200µl
Material Characteristics	polypropylene (PP), premium surface
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile R SAL 10-6
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens and endotoxins
	Sterile R
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months
Packaging	pcs per inner carton: 960
	pcs per carton: 5.760
	pcs per pallet: 172.800
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

Technical Datasheet

Date: 02.01.25

Page 2 of 2

Approved by nerbe plus quality department