

Technical Datasheet

Date: 15.08.25

Page 1 of 2

EN ISO 13485 certified		
Article - No:	06-695-5300	Revision: 7
Remarks		single use only

Specification

Product	Filter tip PP, premium surface, 100-1.250µl, 96 pcs/rack, transparent,, moulded rings, np premium bio, sterile R
Dimensions	length: 98,20 mm (+/- 0,2 mm)
	inner diameter: 0,90 mm (+/- 0,2 mm)
Colour	transparent,
Graduation	moulded rings
Nominal Capacity	100-1.250µl
Material Characteristics	polypropylene (PP), premium surface
	filter: high density polyethylene (PE)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile R SAL 10-6
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens and endotoxins
	Sterile R
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months
Packaging	pcs per inner carton: 576
	pcs per carton: 3.456
	pcs per pallet: 103.680
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

Technical Datasheet

Date: 15.08.25

Page 2 of 2

Approved by nerbe plus quality department