

PŘEHLED TESTŮ RIDA®GENE A KOMPATIBILNÍCH SYSTÉMŮ qPCR

		Real-time PCR přístroje				
		RIDA®CYCLER (R-Biopharm AG)	LightCycler® 480 II (Roche)	ABI 7500 Fast (Applied Biosystems) ¹	CFX96™ (Bio-Rad)	Rotor-Gene Q (Qiagen)
Gastrointestinální infekce						
Viry						
PG1315	RIDA®GENE Viral Stool Panel I		✓			
PG1325	RIDA®GENE Viral Stool Panel II	✓	✓	✓	✓	✓
PG1335	RIDA®GENE Viral Stool Panel III	✓	✓	✓	✓	✓
PG1405	RIDA®GENE Norovirus	✓	✓	✓	✓	✓
PG1415	RIDA®GENE Norovirus I & II	✓	✓	✓	✓	✓
PG1605	RIDA®GENE Sapovirus	✓	✓	✓	✓	✓
PG1605	RIDA®GENE Enterovirus	✓	✓	✓	✓	✓
Bakterie						
PG2205	RIDA®GENE EHEC/EPEC	✓	✓	✓	✓	✓
PG2215	RIDA®GENE EAEC	✓	✓	✓	✓	✓
PG2225	RIDA®GENE ETEC/EIEC	✓*	✓	✓	✓	✓
PG2285	RIDA®GENE E. coli Stool Panel I	✓	✓	✓	✓	✓
PG2305	RIDA®GENE Helicobacter pylori	✓	✓	✓	✓	✓
PG2405	RIDA®GENE Bacterial Stool Panel	✓	✓	✓	✓	✓
PG2415	RIDA®GENE Bacterial Stool Panel I		✓			
Paraziti						
PG1715	RIDA®GENE Parasitic Stool Panel I		✓			
PG1725	RIDA®GENE Parasitic Stool Panel II	✓	✓	✓	✓	✓
PG1745	RIDA®GENE Dientamoeba fragilis	✓	✓	✓	✓	✓
Nozokomiální infekce						
Viry a bakterie						
PG0605	RIDA®GENE MRSA	✓	✓	✓	✓	✓
PG0705	RIDA®GENE Hospital Stool Panel	✓	✓	✓	✓	✓
PG0835	RIDA®GENE Clostridium difficile	✓	✓	✓	✓	✓
Gastrointestinální infekce						
Viry						
PG0505	RIDA®GENE Flu	✓	✓	✓	✓	✓
PG0545	RIDA®GENE Flu & RSV	✓	✓	✓	✓	✓
PG1005	RIDA®GENE Adenovirus	✓	✓	✓	✓	✓
PG4725	RIDA®GENE Pediatric Viral Panel	✓	✓	✓	✓	✓
PG5805	RIDA®GENE Parainfluenza	✓	✓	✓	✓	✓
PG5905	RIDA®GENE RSV & hMPV	✓	✓	✓	✓	✓
PG6805	RIDA®GENE Coronavirus	✓	✓	✓	✓	✓
PG6815	RIDA®GENE SARS-CoV-2	✓	✓	✓	✓	✓
PG6825	RIDA®GENE Flu & SARS-CoV-2	✓	✓	✓	✓	✓

✓ - k dispozici

✓ - brzy k dispozici

* - odchýlný detekční limit cyklů

PŘEHLED TESTŮ RIDA® GENE A KOMPATIBILNÍCH SYSTÉMŮ qPCR

		Real-time PCR přístroje				
		RIDA®CYCLER (R-Biopharm AG)	LightCycler® 480 II (Roche)	ABI 7500 Fast (Applied Biosystems)¹	CFX96™ (Bio-Rad)	Rotor-Gene Q (Qiagen)
Respirační infekce						
Bakterie						
PG2505	RIDA®GENE Bordetella	✓	✓	✓	✓	✓
PG2705	RIDA®GENE CAP Bac	✓	✓	✓	✓	✓
PG8005	RIDA®GENE Legionella	✓	✓	✓	✓	✓
Plísňe						
PG1905	RIDA®GENE Pneumocystis jirovecii		✓	✓	✓	✓
Pohlavně přenosné infekce						
Bakterie						
PG4945	RIDA®GENE STI Mycoplasma Panel	✓	✓	✓	✓	✓
Infekce kůže a měkkých tkání						
Bakterie						
PG0645	RIDA®GENE PVL	✓	✓	✓	✓	✓
Narušená bakteriální rovnováha						
Bakterie						
PG0105	RIDA®GENE Gut Balance		✓	✓	✓	✓
PG0145	RIDA®GENE Akkermansia muciniphila		✓	✓	✓	✓
PG0155	RIDA®GENE Faecalibacterium prausnitzii		✓	✓	✓	✓
Lidská genetika						
PY0205	RIDA®GENE HLA-B27		✓	✓	✓	✓
PY1205	RIDA®GENE Factor II		✓	✓	✓	
PY1210	RIDA®GENE Factor V		✓	✓	✓	
PY4215	RIDA®GENE Lac Intol		✓	✓	✓	

✓ - k dispozici



BioVendor - Laboratorní medicína s.r.o.

Karásek 1767/1, 621 00 Brno

Česká republika

www.biovendor.cz

Kontaktujte nás

Mgr. Eva Dulavová

+420 606 646 629

dulavova@biovendor.cz

