

Liofilchem®

Certificate of Analysis

Page 1 of 1

Product	Batch	Expiration date
Piperacillin - tazobactam T2P 36 µg	050826101	2029.05.07

Ref.

9160 – 9160/1 – 9160/2

Antimicrobial Susceptibility Testing

Tested according to current EUCAST methodology

Control strains	Medium	Inoculum	Incubation	Expected Results Zone range (mm)	Results Zone (mm)
<i>Escherichia coli</i> ATCC® 25922	Mueller Hinton II Agar	0.5 McFarland	35 ± 1°C, ambient 16-18 h	21-27	27
<i>Pseudomonas aeruginosa</i> ATCC® 27853	Mueller Hinton II Agar	0.5 McFarland	35 ± 1°C, ambient 16-18 h	23-29	28
<i>Escherichia coli</i> ATCC® 35218	Mueller Hinton II Agar	0.5 McFarland	35 ± 1°C, ambient 16-18 h	21-27	27
<i>Klebsiella pneumoniae</i> ATCC® 700603	Mueller Hinton II Agar	0.5 McFarland	35 ± 1°C, ambient 16-18 h	14-20	19
<i>Haemophilus influenzae</i> ATCC® 49766	Mueller Hinton Fastidious Agar	0.5 McFarland	35 ± 1°C, 5% CO ₂ 16-18 h	32-40	36
<i>Bacteroides fragilis</i> ATCC® 25285	Fastidious Anaerobe Agar w/ Horse Blood	1 McFarland in broth	36 ± 1°C, anaerobically, 24-72 h	29-35	32
<i>Clostridium perfringens</i> ATCC® 13124	Fastidious Anaerobe Agar w/ Horse Blood	1 McFarland in broth	36 ± 1°C, anaerobically, 24-72 h	29-35	32

Batch Release

Approved

Date 18.05.2026

Signature

Quality Control
(This document has been
established electronically and
is valid without signature)

The results reported were obtained at the time of release.