

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 700 144-0  
F: +49 30 700 144-30  
info@sifin.de  
www.sifin.de

|                    |                                                                                                                                  |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b>     | <b>Anti-Salmonella O:9</b>                                                                                                       |
| Catalogue no.      | TR1307-01                                                                                                                        |
| Package size       | 5 ml                                                                                                                             |
| <b>Lot no.</b>     | <b>6190824</b>                                                                                                                   |
| <b>Expiry date</b> | <b>2029-07</b>                                                                                                                   |
| Intended use       | Anti-Salmonella O:9 records all <i>Salmonella</i> strains, as long as they contain antigen O:9 (group D) in slide agglutination. |
| Specificity        | No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.                                |
| Preservatives      | Sodium azide 0.9 mg/ml                                                                                                           |
| Storage            | 2 °C...8 °C                                                                                                                      |

This document has been established electronically and is valid without signature.  
The batch release was done in accordance with SOP QMV 7.