



LI 1000

Working Volume	With Screw Caps	1,000 μ l
	With Push Caps	1,070 μ l
	With Heat Seal	1,120 μ l
Type of Rack	Polycarbonate Polypropylene	
Type of Rack Lid	Slide Lock Push Lock	
Rack Identifier	Narrow Side Bottom	1D and Human Readable 2D
Storage Temperature	With Screw Caps	-196°C
	With Push Caps	-80°C
	With Heat Seal	-80°C
Standard Side Code Tube	(optional)	1D and Human Readable
Screw Caps	Pre-Capped Stacked on SBS 96 Cap Rack Bulk Ware	
Purity	Production Sterile (ISO 8) Pre-Sterilized (Beta Radiation)	
Outer Diameter	8.65 mm	
Tube Height	Without Cap	44.8 mm
	With Screw Cap	51.2 mm
	With Push Cap	45.9 mm
Rack Height	Without Cap	45.3 mm
	Screw Capped Tubes	51.7 mm
	With Lid	54.3 mm



Technical Drawing



Standard Item Number 2DNC-I10-BL-NS-SLP-L

Order Information

