
SX 300 Mini

Working Volume	With Screw Caps	280 μ	ı
	With Push Caps	280 μ	l
	With Heat Seal	310 µ	i
Type of Rack	Polycarbonate		
Type of Rack Lid	Slide Lock		
Rack Identifier	Narrow Side	1D an	d Human Readable
	Bottom	2D	
Storage Temperature	With Screw Caps	-196°(
	With Push Caps	-80°C	
	With Heat Seal	-80°C	
Standard Side Code Tube	(optional)	Huma	n Readable
Screw Caps	Pre-Capped		
	Stacked on SBS 96 Cap Rack		
	Bulk Ware		
Purity	Production Sterile (ISO 8)		
	Pre-Sterilized (Beta Radiation)		
Outer Diameter	8.65 mm		
Tube Height	Without Cap	17.0 m	nm
	With Screw Cap	18.4 n	nm
	With Push Cap	18.4 n	nm
Rack Height	Without Cap	1	7.5 mm
	Screw Capped Tu	bes 1	8.9 mm
			2.3 mm



Technical Drawing



Standard Item Number 2DNC-X028BL-NS-SLC-L

Order Information

