

# SAFE® CAP 8 CHANNEL DD

## Automated 8-Channel Capper/Decapper

- Fully-automated capper/decapper as desktop device
- Reliably decaps and caps eight tubes in parallel
- Stable stand for storage of device and contamination-free parking of screw caps
- Eight individual motors for precise application of 7 Ncm torque for tube closure
- Also integrable in automated pipetting solutions

### Item Number


SAFE® 96 XT:

SAFE® 96 IT:

CDC-8CH-XT-DD

CDC-8CH-IT-DD

### Order Information



SAFE® CAP 8 CHANNEL DD