

SAFE[®] ACCESS TUBE PICKER

Semi-automated Tube Picker

- Fast and reliable tube picking from 96 SBS racks
- Works with frozen and thawed tubes
- Generate picklists externally in other software applications, or create a template using the own software
- Can be connected to a 2D barcode scanner for tube confirmation and for generating picklists directly from the tube barcodes
- 1D scanner for rack barcode included as standard
- Tubes can be picked from a single rack or a picklist can be set up to select tubes across multiple racks
- Picklist import file formats are Excel, Text, JSON and XML
- Operating system on host PC = Windows 7, 8, 10
- Dimensions: 300 mm x 300 mm x 450 mm

Item Number

TPI96

Order Information

