

CERTIFICATE OF ANALYSIS

		MIC_COA_SI390900#SCR23304.pdf				
Product Code	Name	Quantity manufactured	d			
SI 390.900	URO-QUICK SCREENING KI	27				
Lot Number	Manufacturing Date	Expiry Date	CAP color	Printed Band		
SCR23304	07/09/2023	05/11/2024	Green	White		

Eugonic Broth

Lot Number SCR23304		Physical parameters (Volume, transparency, colour)	Decontamination at 37°C after 72 hours	Growth inter batch	Growth intra batch
Expiry Date 05/11/2024	Expected	Compliant	No growth	CV%<=35%	CV%<=35%
CAP color Green Printed Band White	Found	ОК	ОК	17.8%	< 5%

Check of the presenze of the magnetic stirrer in each vial contained in the kit: YES

Storage Conditions: min 4°C - Max 30°C

Compliance: YES

Controlled By	Cirio Alberto	Date of Control	13/09/2023
Approved By Manufacturing Manager	Pividori Denis	Date of Approval	13/09/2023

We certify that the lot in object has been manufactured and tested according with ALIFAX manufacturing and quality control procedures and specifications. We confirm that the above referenced product conformed to the specifications and then it is approved.

Approved by Quality Assurance Falcone Francesco

Date of Release 14/09/2023

This is a computer-generated document. No signature is required.

Registered Office

Via Petrarca, 2/1 - 35020 Polverara (PD) Italy Tel (+39) 049 0992000 - Fax (+39) 049 5855434 Mail PEC: info.alifax@pec.alifax.com info@alifax.com - www.alifax.com



Manufacturing Site

Via Merano, 30 - 33045 Nimis (UD) Italy Tel (+39) 0432 547454 - Fax (+39) 0432 547378 Mail PEC: info.udine@pec.alifax.com info.udine@alifax.com