

CERTIFICATE OF ANALYSIS

			MIC_COA_SI390900#SCR21018.pdf			
Product Code	Name	Quantity manufactured				
SI 390.900	URO-QUICK SCREENING KIT	27				
Lot Number	Manufacturing Date	Expiry Date	CAP color	Printed Band		
SCR21018	25/01/2021	26/03/2022	Green	White		

Eugonic Broth

Lot Number	SCR21018		Physical parameters (Volume, transparency, colour)	Decontamination at 37°C after 72 hours	Growth inter batch	Growth intra batch
Expiry Date	26/03/2022	Expected	Compliant	No growth	CV%<=35%	CV%<=35%
CAP color	Green	Found	ОК	ОК	21.9%	23.6%
Printed Band	White					

Check of the presenze of the magnetic stirrer in each vial contained in the kit: YES

Storage Conditions: min 4°C - Max 30°C

YES

Compliance:

Controlled By	Cirio Alberto	Date of Control	29/01/2021
Approved By Manufacturing Manager	Pividori Denis	Date of Approval	29/01/2021

We certify that the lot in object has been manufactured and tested according with ALIFAX manufacturing and quality control procedures and specifications. These requirements have been established according to the European Directive 98/79/EC. We confirm that the above referenced product conformed to the specifications and then it is released.

Released By QualityFacchin FabioDate of Release29/01/2021

This is a computer-generated document. No signature is required.

Registered Office

Via Petrarca, 2/1 - 35020 Polverara (PD) Italy Tel (+39) 049 0992000 - Fax (+39) 049 5855434 Mail PEC: info.alifax@pec.alifax.com info@alifax.com - www.alifax.com



Manufacturing Site

Via Merano, 30 - 33045 Nimis (UD) Italy Tel (+39) 0432 547454 - Fax (+39) 0432 547378 Mail PEC: info.udine@pec.alifax.com info.udine@alifax.com