


Certificate of Quality Control

| | |
|---------------------------------|--|
| Manufacturer | SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany |
| Product | Test Serum Anti-Yersinia enterocolitica O 5 |
| Catalogue No. (Package size) | TS 1704 (1 ml) |
| Lot No. | 158 11 12 |
| Expiry date | 2017-11 |
| Intended use | Records all Yersinia enterocolitica strains in the slide agglutination as long as they contain the O 5 antigen |
| Specificity | No agglutination with heterologous Yersinia enterocolitica O antigens and selected representatives of other Enterobacteriaceae. |
| Preservatives | Sodium azide 0.9 mg / ml |
| Storage | 2...8 °C |

Berlin, 11. 12. 12


Head of Quality Control