

### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Yersinia enterocolitica-O9 Test Antigen for the Micro Widal Reaction</b>
Catalogue No. (Package size)	TS 1723 (3 ml)
Lot No.	<b>71 06 13</b>
Expiry date	<b>2018-06</b>
Intended use	In the micro Widal test, the test antigen reacts with the homologous <i>Yersinia enterocolitica</i> O9 antibodies in the serum to be tested. The specific bonding is indicated by agglutination.
Specificity	There must be no reaction with heterologous antibodies anti- <i>Yersinia</i> O3, O5 and O6 or with physiological saline solution.
Preservatives	Thiomersal 0.2 mg/ml
Storage	2...8 °C

Berlin,

10.06.13

  
Head of Quality Control