

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Coli I</b>
Catalogue no.	TS 2111
Package size	1 ml
<b>Lot no.</b>	<b>2670317</b>
<b>Expiry date</b>	<b>2022-03</b>
Intended use	Anti-Coli I records all <i>E. coli</i> strains, as long as they contain the antigens O 26:(K 60), O 44:(K 74), O 114:(K 90), O 125:(K 70), O 142:(K 86) and O 158:(K-).
Specificity	No agglutination with <i>E. coli</i> strains which are recorded by Anti-Coli II and Anti-Coli III.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2 °C...8 °C

Berlin,  
place, date

6.07.2017

  
signature Head of Quality Control