

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Anti-Coli I
Catalogue no.	TS 2111
Package size	1 ml
Lot no.	0071119
Expiry date	2024-11
Intended use	Anti-Coli I records all <i>E. coli</i> strains, as long as they contain the antigens O 26:(K 60), O 44:(K 74), O 114:(K 90), O 125:(K 70), O 142:(K 86) and O 158:(K-).
Specificity	No agglutination with <i>E. coli</i> strains which are recorded by Anti-Coli II and Anti-Coli III.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2 °C...8 °C

Berlin, 12.12.19
place, date

F. Leuch
signature Head of Quality Control