

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella H:z
Catalogue no.	TR1424
Package size	1 ml
Lot no.	8541121
Expiry date	2025-04
Intended use	Anti-Salmonella H:z records all <i>Salmonella</i> strains in Subspecies I, II, III, as long as they contain antigen H:z in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 06.12.21
place, date

F. Lendy
signature Head of Quality Control