

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

| | |
|--------------------|--|
| Product | Anti-Yersinia enterocolitica O 3 |
| Catalogue no. | TS1701 |
| Package size | 1 ml |
| Lot no. | 2130322 |
| Expiry date | 2027-02 |
| Intended use | Records all <i>Yersinia enterocolitica</i> strains as long as they contain antigen O 3 in slide agglutination. |
| Specificity | No agglutination with heterologous <i>Yersinia enterocolitica</i> O antigens and selected representatives of other <i>Enterobacteriaceae</i> . |
| Preservatives | Sodium azide 0.9 mg/ml |
| Storage | 2 °C...8 °C |

Berlin,

19.05.22

place, date



signature Head of Quality Control