

Liofilchem®

Certificate of Analysis

Page 1 of 1

| Product | Batch | Expiration date |
|-----------------------|-----------|-----------------|
| Fosfomycin FOS 200 µg | 031324055 | 2027.03.13 |

Ref.

9109 – 9109/1 – 9109/2

Antimicrobial Susceptibility Testing

Tested according to current CLSI and/or EUCAST methodology

| Control strains | Medium | Inoculum | Incubation | Expected Results Zone range (mm) | Results Zone (mm) |
|---|---------------------------|---------------|------------------------------|-------------------------------------|----------------------|
| <i>Escherichia coli</i> ATCC® 25922 | Mueller Hinton II Agar | 0.5 McFarland | 35 ± 1°C, ambient 16-18 h | 22-30 | 27 |
| <i>Escherichia coli</i> ATCC® 25922 | Mueller Hinton II Agar | 0.5 McFarland | 35 ± 1°C, ambient 16-18 h | 26-34* | 27 |
| <i>Staphylococcus aureus</i> ATCC® 25923 | Mueller Hinton II Agar | 0.5 McFarland | 35 ± 1°C, ambient 16-18 h | 25-33 | 28 |

*Established and validated by EUCAST

Batch Release

Approved

Date

20.03.2024

Signature

Quality Control

(This document has been
established electronically and
is valid without signature)

The results reported were obtained at the time of release.