


Certificate of Analysis

Product specifications

Product: CF StripAssay®
Catalog Number: 4-410
 4-410-TRIAL
Lot: 10-EI-26-036
Date of Manufacture: 2026-03-26
Expiry Date: 2027-03
Storage temperature: 2-8°C

Quality Control Testing

Parameter	Tested	Corresponds
Specificity/Sensitivity		
Detectability of each allele/genotype	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Lack of cross-hybridization	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Reproducibility		
Conformity of results in multiple runs	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Lack of nucleic acids contamination	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Physical parameters (color, pH,...)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

This product complies with the General Safety and Performance Requirements of the Regulation 2017/746/EU on *in vitro* diagnostic medical devices and is marked with the  0123 - symbol.



Manager QC/QA

26.03.2026
Date