


## Certificate of Analysis

### Product specifications

**Product:** CF StripAssay<sup>®</sup>  
**Catalog Number:**  4-410  
 4-410-TRIAL  
**Lot:** 20-ZW-26-016  
**Date of Manufacture:** 2026-04-09  
**Expiry Date:** 2027-05  
**Storage temperature:** 2-8°C

### Quality Control Testing

Parameter	Tested	Corresponds
Specificity/Sensitivity		
Detectability of each allele/genotype	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Lack of cross-hybridization	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Reproducibility		
Conformity of results in multiple runs	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Lack of nucleic acids contamination	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Physical parameters (color, pH,...)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

This product complies with the General Safety and Performance Requirements of the Regulation 2017/746/EU on *in vitro* diagnostic medical devices and is marked with the  0123 - symbol.



Manager QC/QA

09.04.2026  
Date