

Solstice

CERTIFICATE OF ANALYSIS

Honeywell Specialty Chemicals Seelze GmbH D-30926

Wunstorfer Str. 40

Seelze, 30926

Telefon: +49 5137 999-0

lab.honeywell.com

Brand:

Product:

Lot No:

Production date:

Expiry Date:

Honeywell Riedel-de Haën

34863

P266A

23.Sep.2025

06.Mar.2029

2-Propanol CHROMASOLV, 99.9%

Parameter	Specification	Units	Result	Units
IR spectrum	complying		complying	
assay (GC)	min. 99.85	%	>99.99	%
non-volatile matter	max. 0.0003	%	<0.0003	%
water (Karl Fischer)	max. 0.05	%	0.006	%
free acid (as C2H5COOH)	max. 0.001	%	<0.001	%
absorbance at 400 nm	max. 0.01		<0.01	
absorbance at 260 nm	max. 0.01		<0.01	
absorbance at 230 nm	max. 0.10		<0.10	
absorbance at 220 nm	max. 0.20		0.20	
absorbance at 210 nm	max. 0.70		<0.70	
absorbance at 205 nm	max. 1.0		<1.0	
APHA	max. 10		<10	

Solstice

filtered via 0,2µm pore size filter

QC Release Date: 06.Oct.2025

The minimum shelf life is based on the current knowledge and holds only for proper storage conditions in the originally closed flasks/ packages.

We herewith confirm that the delivery is effected according to the technical delivery conditions agreed.

Particular properties of the products or the suitability for a particular area of application are not assured.

We guarantee a proper quality within our General Conditions of Sales.

Honeywell
Specialty Chemicals

This document was produced electronically and is valid without a signature