



CERTIFICATE OF ANALYSIS

Honeywell Specialty Chemicals Seelze GmbH
Wunstorferstrasse 40
Seelze, 30926
Telefon: +49 5137 999-0
www.lab-honeywell.com

Brand: Honeywell Riedel-de Haën
Product: 34863
Lot No: L007B
Production date: 07.Jan.2021
Expiry Date: 20.Jun.2024

2-Propanol CHROMASOLV, 99.9%

Parameter	Specification	Units	Result	Units
IR spectrum	complying		complying	
assay (GC)	min. 99.85	%	>99.99	%
non-volatile matter	max. 0.0003	%	<0.0003	%
water (Karl Fischer)	max. 0.05	%	0.002	%
free acid (as C ₂ H ₅ COOH)	max. 0.001	%	<0.001	%
absorbance at 400 nm	max. 0.01		<0.01	
absorbance at 260 nm	max. 0.01		<0.01	
absorbance at 230 nm	max. 0.10		<0.10	
absorbance at 220 nm	max. 0.20		<0.20	
absorbance at 210 nm	max. 0.70		<0.70	
absorbance at 205 nm	max. 1.0		<1.0	
APHA	max. 10		<10	



filtered via 0,2µm pore size filter

QC Release Date: 28.Jan.2021

The minimum shelf life is based on the current knowledge and holds only for proper storage conditions in the originally closed flasks/ packages.

We herewith confirm that the delivery is effected according to the technical delivery conditions agreed.

Particular properties of the products or the suitability for a particular area of application are not assured.

We guarantee a proper quality within our General Conditions of Sales.

Honeywell

Quality Management

This document was produced electronically and is valid without a signature