

**ELI.H.A *Aspergillus* (102 tests)**

Référence / Catalog number : 44602

Lot de conditionnement / Packaging lot : 8939

Date d'expiration / Expiration date : 2016-12-31

Réactifs / Reagents	Lot de fabrication / Production lot	Date d'expiration / Expiration date
Flacon R1 (2,2 mL) / R1 (2.2 mL) vial	8939	2016-12-31
Flacon R2 (1 mL) / R2 (1 mL) vial	8939	2016-12-31
Flacon CONTROL + (0,2 mL) / CONTROL + (0.2 mL) vial	8939	2016-12-31
Flacon CONTROL - (0,2 mL) / CONTROL - (0.2 mL) vial	8939	2016-12-31
Flacon R3 (2 mL) / R3 (2 mL) vial	8939	2016-12-31
Flacon BUF (55 mL) / BUF (55 mL) vial	8939	2016-12-31

Certificat de contrôle établi suivant les procédures de contrôles internes  
*Control certificate established according to the internal control procedures*

	Résultats / Results	Résultats théoriques / Theoretical results
<b>R1</b> Flacon d'Hématies sensibilisées <i>Vial of Sensitized red blood cells</i>		
Conditionnement primaire / Primary packaging	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
Etiquetage / Labelling	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
Aspect du milieu / Medium appearance	Conforme / <i>Conform to specifications</i>	Suspension d'hématies marron <i>Brown suspension of red blood cells</i>
Performances / Performances	Conforme / <i>Conform to specifications</i>	Sérums humains frais / Fresh human sera Sérums positifs / Positive sera Conformité / spécifications <i>Conformity / specifications</i>
<b>R2</b> Flacon d'Hématies non sensibilisées <i>Vial of Unsensitized red blood cells</i>		
Conditionnement primaire / Primary packaging	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
Etiquetage / Labelling	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
Aspect du milieu / Medium appearance	Conforme / <i>Conform to specifications</i>	Suspension d'hématies marron <i>Brown suspension of red blood cells</i>
Performances / Performances	Conforme / <i>Conform to specifications</i>	Sérums humains frais / Fresh human sera Sérums positifs / Positive sera Conformité / spécifications <i>Conformity / specifications</i>

	Résultats / Results	Résultats théoriques / Theoretical results
<b>CONTROL +</b> Flacon de Contrôle positif <i>Vial of Positive control</i>		
Conditionnement primaire / <i>Primary packaging</i>	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
Etiquetage / <i>Labelling</i>	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
Aspect du milieu / <i>Medium appearance</i>	Conforme / <i>Conform to specifications</i>	Liquide jaunâtre translucide <i>Translucent yellowish liquid</i>
Titre / <i>Titer</i>	1 / 320 Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
<b>CONTROL -</b> Flacon de Contrôle négatif <i>Vial of Negative control</i>		
Conditionnement primaire / <i>Primary packaging</i>	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
Etiquetage / <i>Labelling</i>	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
Aspect du milieu / <i>Medium appearance</i>	Conforme / <i>Conform to specifications</i>	Liquide jaunâtre translucide <i>Translucent yellowish liquid</i>
<b>R3</b> Flacon d'adsorbant <i>Vial of Adsorbent</i>		
Conditionnement primaire / <i>Primary packaging</i>	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
Etiquetage / <i>Labelling</i>	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
Aspect du milieu / <i>Medium appearance</i>	Conforme / <i>Conform to specifications</i>	Suspension d'hématies marron <i>Brown suspension of red blood cells</i>
<b>BUF</b> Flacon de tampon phosphate <i>Vial of Phosphate buffer</i>		
Conditionnement primaire / <i>Primary packaging</i>	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
Etiquetage / <i>Labelling</i>	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>
Aspect du milieu / <i>Medium appearance</i>	Conforme / <i>Conform to specifications</i>	Liquide incolore translucide <i>Translucent colorless liquid</i>
<b>Contrôle du conditionnement / <i>Conformity after packaging</i></b>	Conforme / <i>Conform to specifications</i>	Conformité / spécifications <i>Conformity / specifications</i>

Référence / *Catalog number* : 44602  
 Lot de conditionnement / *Packaging lot* : 8939

Lot accepté le / *Accepted on* : 21/07/2015

Par / *By* : L.ROBIN  
 Responsable Contrôle Qualité  
 Quality Manager

