

Liofilchem®

Certificate of Analysis

Page 1 of 1

| Product | Batch | Production date | Expiration date |
|---------------------------------|-----------|-----------------|-----------------|
| Teicoplanin TEC 0.016-256 µg/mL | 032017060 | 20.03.2017 | 2020.03.19 |

Ref.

92012 – 920120 – 920121

Antimicrobial Susceptibility Testing

Tested according to current CLSI methodology

| Control strains | Medium | Inoculum | Incubation | Expected Results MIC range (mg/L) | Results MIC(mg/L) |
|---|---------------------------|----------------------------|------------------------------|--------------------------------------|----------------------|
| <i>Staphylococcus aureus</i> ATCC® 29213 | Mueller Hinton II Agar | 0.5 McFarland in saline | 35 ± 2°C, ambient 16-20 h | 0.25–1 | 0.75 |
| <i>Enterococcus faecalis</i> ATCC® 29212 | Mueller Hinton II Agar | 0.5 McFarland in saline | 35 ± 2°C, ambient 16-20 h | 0.25–1 | 0.5 |

Batch Release

Approved

Date

26.04.2017

Signature

Quality Control
(D. Vitagliano)



The results reported were obtained at the time of release.