



www.pro-lab.com

CERTIFICATE OF ANALYSIS

Product Name: Prolex™ <i>E. coli</i> O157 Latex Test Reagent Kit		Lot Number: C11577
Product Code: PL.071B		Expiration Date: 2015 09
Product Code	Quantity	Components
PL.073B	1 x 6.2 ml	<i>E. coli</i> O157 Latex Reagent
PL.075B	1 x 3.0 ml	<i>E. coli</i> O157 Positive Control
PL.076B	1 x 3.0 ml	<i>E. coli</i> O157 Latex Reagent Negative Control
	150	Mixing Sticks
PL.092-20	2 x 20 cards	Test Cards
Additional Information		
Storage Conditions:	2-8°C	
Precautions:	Refer to Instructions for Use	
Quality Control Testing		
	Criteria	Result
	Strong agglutination of Latex Reagent when tested against <i>E. coli</i> serotype O157 isolates.	PASS
	No agglutination reaction of Latex Reagent when tested against non O157 serotypes <i>E. coli</i> , and <i>E. hermanii</i> isolates.	PASS
	No agglutination reaction of Negative Control when tested against <i>E. coli</i> and <i>E. hermanii</i> isolates.	PASS
	Strong agglutination of the <i>E. coli</i> O157 Positive Control when tested against Latex Reagent.	PASS
	No agglutination of the <i>E. coli</i> O157 Positive Control when tested against Negative Control.	PASS
Release		
	2014 04 16	<i>Shilpa Tejpal</i>
	Date	Quality Control

Canada ☎ 800 268 2341 Fax 905 731 0206
 20 Mural Street, Unit #4, Richmond Hill, Ontario,
 L4B 1K3
 support@pro-lab.com

U.S.A. ☎ 800 522 7740 Fax 800 332 0450
 21 Cypress Blvd., Suite 1070, Round Rock,
 Texas 78665-1034
 ussupport@pro-lab.com

U.K. ☎ 0151 353 1613 Fax 0151 353 1614
 3 Bassendale Road, Bromborough, Wirral,
 Merseyside CH62 3QL
 uksupport@pro-lab.com