



CERTIFICATE OF ANALYSIS

Product Name: Proflow Legionella Urinary Antigen	Lot Number: LN-084
Product Code: PL.3007	Expiration Date: 2022 02 28

Additional Information

Test conditions: As defined in the IFU. Storage conditions: 2-30°C The product is intended for in vitro use only.

Quality Control Testing

Functionality test results: visual interpretation at 10 minutes

- 1. Sensitivity test.
Dilutions from standard: approved.
- 2. Specificity test.
Materials prepared in product buffer: approved.

Quality Control result

Approved.

Release

2020/02/24	<i>mawr</i>
Date	Quality Control

Canada ☎ 800 268 2341 Fax 905 731 0206
20 Mural Street, Unit #4, Richmond Hill, Ontario,
L4B 1K3
support@pro-lab.com

U.S.A. ☎ 800 522 7740 Fax 800 332 0450
21 Cypress Blvd., Suite 1070, Round Rock,
Texas 78665-1034
ussupport@pro-lab.com

U.K. ☎ 0151 353 1613 Fax 0151 353 1614
3 Bassendale Road, Bromborough, Wirral,
Merseyside CH62 3QL
uksupport@pro-lab.com