



CERTIFICATE OF ANALYSIS

ESR_COA_SI305#1704A.pdf

Product Code SI 305.300-A	Name LATEX CONTROLS KIT - provetta Greiner30 TESTS	Quantity manufactured 178	Type of Tubes Greiner
-------------------------------------	--	-------------------------------------	---------------------------------

Lot Number 1704A	Manufacturing Date 01/07/2016	Expiry Date 28/12/2016		
----------------------------	---	----------------------------------	--	--

Standard Values			Measured Values				Internal	Volume	
Transmittance	ESR (mm/h) (min-Val-max)		ESR (mm/h) (min-Val-max)		Instrument	Serial Number	OK	Aggregation	checked
Level 2	1534	6 - 11	9 9 9	TEST 1 BCL	2173	<input checked="" type="checkbox"/>	NO	YES	
Level 3	980	15 - 22	19 19 20	TEST1 THL	PRDX	<input checked="" type="checkbox"/>	NO	YES	
Level 4	637	56 - 74	65 66 66	ROLLER 20 PN	133	<input checked="" type="checkbox"/>	NO	YES	

Date of Control **01/07/2016**

Date of Approval **01/07/2016**

Controlled By **Pividori Denis**

Approved By **Ferino Rudy**

Registered Office

Via Petrarca, 2/1 - 35020 Polverara (PD) Italy
Tel (+39) 049 0992000 - Fax (+39) 049 5855434
Mail PEC: info.alifax@pec.alifax.com
info@alifax.com - www.alifax.com



Manufacturing Site

Via Merano, 30 - 33045 Nimis (UD) Italy
Tel (+39) 0432 547454 - Fax (+39) 0432 547378
Mail PEC: info.udine@pec.alifax.com
info.udine@alifax.com