



CERTIFICATE OF ANALYSIS

ESR_COA_SI305#2164B.pdf

Product Code SI 305.300-A	Name LATEX CONTROLS	Quantity manufactured 176	Type of Tubes Greiner
-------------------------------------	-------------------------------	-------------------------------------	---------------------------------

Lot Number 2164B	Manufacturing Date 29/01/2021	Expiry Date 28/07/2021		
----------------------------	---	----------------------------------	--	--

Standard Values			Measured Values				Internal Aggregation	Volume checked	
Transmittance	ESR (mm/h) (min-Val-max)		ESR (mm/h) (min-Val-max)			Instrument	Serial Number	OK	
Level 2	1534	6 - 11	9	9	9	TEST 1 BCL	2173	<input checked="" type="checkbox"/>	NO YES
Level 3	980	15 - 22	19	19	20	TEST1 THL	T13290	<input checked="" type="checkbox"/>	NO YES
Level 4	637	56 - 74	65	66	68	ROLLER 20 PN	133	<input checked="" type="checkbox"/>	NO YES

Controlled By

Cirio Alberto

Date of Control

29/01/2021

Approved By Manufacturing Manager

Pividori Denis

Date of Approval

29/01/2021

We certify that the lot in object has been manufactured and tested according with ALIFAX manufacturing and quality control procedures and specifications. These requirements have been established according to the European Directive 98/79/EC. We confirm that the above referenced product conformed to the specifications and then it is released.

Released By Quality

Facchin Fabio

Date of Release

29/01/2021

This is a computer-generated document. No signature is required.

Registered Office

Via Petrarca, 2/1 - 35020 Polverara (PD) Italy
Tel (+39) 049 0992000 - Fax (+39) 049 5855434
Mail PEC: info.alifax@pec.alifax.com
info@alifax.com - www.alifax.com



Manufacturing Site

Via Merano, 30 - 33045 Nimis (UD) Italy
Tel (+39) 0432 547454 - Fax (+39) 0432 547378
Mail PEC: info.udine@pec.alifax.com
info.udine@alifax.com