

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella A-67 + Vi, omnivalent
Catalogue no.	TR 1101
Package size	1 ml
Lot no.	3450419
Expiry date	2020-10
Intended use	Anti-Salmonella A-67 + Vi, omnivalent, records all <i>Salmonella</i> serovars of the groups A-67 and the Vi-antigen in slide agglutination.
Specificity	No agglutination with <i>Shigella sonnei</i>
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 16.05.19
place, date

F. Lechner
signature Head of Quality Control