

## Qualitätskontrollzertifikat

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin
Product	<b>Test Reagens Anti-Salmonella OME</b> (O-Group-Pool: X, Y, Z, 51, 52, 53)
Catalogue No. (Package size)	TR 1172 (1 ml)
Lot No.	<b>62 05 13</b>
Expiry date	<b>2015-05</b>
Intended use	Test Reagent Anti-Salmonella OME records all <i>Salmonella</i> strains belonging to the O-Groups X, Y, Z, 51, 52, 53.
Specificity	No agglutination with <i>Salmonella</i> strains which are recorded by Test Reagents Anti-Salmonella OMA (O-Group-Pool: A, B, D, E, L), Anti-Salmonella OMB (O-Group-Pool: C, F, G, H), Anti-Salmonella OMC (O-Group-Pool: I, J, K, M, N, O, P), Anti-Salmonella OMD (O-Group-Pool: Q, R, S, T, U, V, W), Anti-Salmonella OMF (O-Group-Pool: 54-59) and Anti-Salmonella OMG (O-Group-Pool: 60-63, 65-67).
Preservatives	Sodium azide 0,9 mg/ml
Storage	2...8 °C

Berlin,

31.05.13

J. Ludwig  
Head of Quality Control