

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 700 144-0  
F: +49 30 700 144-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella OME</b> O-Group-Pool X, Y, Z, 51, 52, 53
Catalogue no.	TR1172
Package size	1 ml
<b>Lot no.</b>	<b>2870320</b>
<b>Expiry date</b>	<b>2022-02</b>
Intended use	Anti-Salmonella OME records all <i>Salmonella</i> strains belonging to the O-Groups X, Y, Z, 51, 52, 53 in slide agglutination.
Specificity	No agglutination with <i>Salmonella</i> strains which are recorded by Anti-Salmonella OMA (O-Group-Pool: A, B, D, E, L), Anti-Salmonella OMB (O-Group-Pool: C, F, G, H), Anti-Salmonella OMC (O-Group-Pool: I, J, K, M, N, O, P), Anti-Salmonella OMD (O-Group-Pool: Q, R, S, T, U, V, W), Anti-Salmonella OMF (O-Group-Pool: 54-59) and Anti-Salmonella OMG (O-Group-Pool: 60-63, 65-67).
Preservatives	Sodium azide 0,9 mg/ml
Storage	2 °C...8 °C

Berlin, 04.05.2020  
place, date

J. Leuders  
signature Head of Quality Control