

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 700 144-0  
F: +49 30 700 144-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella Group B</b>
Catalogue no.	TR1201
Package size	1 ml
<b>9001019</b>	<b>6590821</b>
<b>Expiry date</b>	<b>2024-07</b>
Intended use	Anti-Salmonella Group B records all <i>Salmonella</i> strains of group B in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,

21.10.20

place, date

*F. Herder*

signature Head of Quality Control