

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella Group C
Catalogue no.	TR 1202
Package size	1 ml
Lot no.	0861217
Expiry date	2020-12
Intended use	Anti-Salmonella Group C records all <i>Salmonella</i> strains of group C (C ₁ ; C ₂ -C ₃) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2°C... 8 °C

Berlin,
place, date

20.01.2018


signature Head of Quality Control