

### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Enteroclon Anti-Salmonella Gr. E</b>
Catalogue No. (Package size)	TR 1204 (1 ml)
Lot No.	<b>127 08 12</b>
Expiry date	<b>2014-08</b>
Intended use	Enteroclon Anti-Salmonella Gr. E records all Salmonella strains of group E in slide agglutination.
Specificity	No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin,

06.08.12

  
\_\_\_\_\_  
Head of Quality Control