

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485


sifin diagnostics gmbh

Berliner Allee 317-321  
 13088 Berlin, Germany

T: +49 30 927 030-0  
 F: +49 30 927 030-30  
 info@sifin.de  
 www.sifin.de

<b>Product</b>	<b>Anti-Salmonella O:6<sub>1</sub></b>
Catalogue no.	TR 1304
Package size	1 ml
<b>Lot no.</b>	<b>73 04 15</b>
<b>Expiry date</b>	<b>2019-04</b>
Intended use	Anti-Salmonella O:6 <sub>1</sub> records all <i>Salmonella</i> strains of the group C <sub>2</sub> -C <sub>3</sub> , as long as they contain antigen O:6 (= O:6 <sub>1</sub> ) in slide agglutination.
Specificity	No agglutination with the serovars of group C <sub>2</sub> -C <sub>3</sub> that contain only the O:8 antigen or the antigen combination O:8,20. The O:6 antigen of the group C <sub>1</sub> (= O: 6 <sub>2</sub> ) is not detected.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 29.04.15  
 \_\_\_\_\_  
 place, date

  
 \_\_\_\_\_  
 signature Head of Quality Control