

## Certificate of Quality Control

sifin diagnostics gmbh

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella O:7</b>
Catalogue no.	TR 1305
Package size	1 ml
<b>Lot no.</b>	<b>80 04 15</b>
<b>Expiry date</b>	<b>2018-04</b>
Intended use	Anti-Salmonella O:7 records all <i>Salmonella</i> strains, as long as they contain antigen O:7 in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C... 8 °C

Berlin,  
place, date

07.05.15

signature Head of Quality Control

