


## Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin
Product	<b>Enteroclon Anti-Salmonella O 9</b>
Catalogue No. (Package size)	TR 1307 (1 ml)
Lot No.	<b>84 09 09</b>
Expiry date	<b>2014-09</b>
Intended use	Enteroclon Anti-Salmonella O 9 records all Salmonella strains, as long as they contain antigen O 9 (group D) in slide agglutination.
Specificity	No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 27.10.08

  
\_\_\_\_\_  
Head of Quality Control