

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin
Product	Enteroclon Anti-Salmonella O 10
Catalogue No. (Package size)	TR 1308 (1 ml)
Lot No.	108 10 09
Expiry date	2014-11
Intended use	Enteroclon Anti-Salmonella O 10 records all Salmonella strains of groups E ₁ and E ₄ , as long as they contain antigen O 10 in slide agglutination.
Specificity	No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 11.12.09

i.V. Ringes
Head of Quality Control