

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

<b>Product</b>	<b>Anti-Salmonella O:10</b>
Catalogue no.	TR1308
Package size	1 ml
<b>Lot no.</b>	<b>0190121</b>
<b>Expiry date</b>	<b>2025-12</b>
Intended use	Anti-Salmonella O:10 records all <i>Salmonella</i> strains of groups E <sub>1</sub> and E <sub>4</sub> , as long as they contain antigen O:10 in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,  
place, date

17.02.21

signature Head of Quality Control

*J. Munding*