


Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella O 19
Catalogue No. (Package size)	TR 1311 (1 ml)
Lot No.	28 02 13
Expiry date	2018-02
Intended use	Enteroclon Anti-Salmonella O 19 records all Salmonella serovars of group E ₄ (O:1,3,19).
Specificity	Enteroclon Anti-Salmonella O 19 cross-reacts slightly with several strains of group B with O 27. However, since the Enteroclon is used only once the Salmonella strain's affiliation to group E has been confirmed, this cross-reaction is irrelevant to the serotyping procedure. No other agglutinations with heterologous Salmonella O antigens and Salmonella H antigens are known.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin,

27.03.13


Head of Quality Control