

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella O 20
Catalogue No. (Package size)	TR 1312 (1 ml)
Lot No.	167 11 12
Expiry date	2016-11
Intended use	Enteroclon Anti-Salmonella O 20 records all Salmonella strains of group C2-C3, as long as they contain antigen O 20.
Specificity	No agglutinations with heterologous Salmonella O antigens and Salmonella H antigens.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 20.12.12


Head of Quality Control