

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

Product	Anti-Salmonella O:20
Catalogue no.	TR 1312
Package size	1 ml
Lot no.	6600718
Expiry date	2022-07
Intended use	Anti-Salmonella O:20 records all <i>Salmonella</i> strains of group C ₂ -C ₃ , as long as they contain antigen O:20 in slide agglutination
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 02.08.2018
 place, date

F. Lueders
 signature Head of Quality Control