

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Enteroclon Anti-Salmonella O 27**

Catalogue No.
(Package size) TR 1313 (1 ml)

Lot No. **111 08 13**

Expiry date **2018-08**

Intended use Enteroclon Anti-Salmonella O 27 records all Salmonella serovars of
group B and D₃ as long as they contain the O 27 antigen.

Specificity No agglutination with heterologous Salmonella O antigens and
Salmonella H antigens.

Preservatives Sodium azide 0.9 mg / ml

Storage 2...8 °C

Berlin,

12.09.13


Head of Quality Control